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

The manufacturer and vendor of this fiberoptic light cable does not accept any liability for direct or consequential damages caused by the improper use or handling, in particular as a result of the failure to follow the instructions below or by improper processing and maintenance.

	The high light concentration at the end of the fiberoptic light cable generates heat. Never place the end of the fiberoptic light cable onto the cover or the skin of the patient, because the light intensity may cause burns to the patient's skin or set the cover on fire.
	Fiberoptic light cables are supplied in non-sterile condition. Therefore, thorough cleaning, disinfection and/or sterilisation are required before the first use as well as before any subsequent utilisation.
	Fiberoptic light cables must not be cleaned in the ultrasonic bath.
	Damage of the coating causes damage to the fiberoptic light cable! Strong kinking causes glass-fiber-breakage and destroys the fiberoptic light cable.
	Damages of the fiberoptic light cables arising as a result of improper handling are not covered by the warranty. It is essential to avoid hitting in particular the adapters and bending the light cables, because this may damage the coating and the optical fibres.
	<p>The cross-sections of the cables must match the diameter of the light input of the endoscope to prevent damages at the light output of the cable caused by reflection.</p> <p>3,5mm - thin endoscopes, such as e.g. for arthroscopies, urological imaging, sinuscopies, etc.</p> <p>4,9mm - thick endoscopes, such as e.g. for laparoscopies, etc.</p> 

Handling and Application

Fiberoptic light cables are used to transmit light energy from a cold light source to an endoscope.

Guideline conformance

This medical device bears the CE mark pursuant to the Medical Device Directive (MDD) 93/42/EEC. Furthermore, every fiberoptic light cable has a separate serial number.

Combination

The RfQ Light Cables are designed for use with halogen or xenon / LED light sources. It is important to note that glued fiber optic cables may only be used in halogen light sources. The following combinations are used.

RfQ-Fiber Light Cable	Connection Light Input	Connection Light Output
00-2xxxx-xx	ACM/Aesculap/Dyonics/ Heine/HSW/Medicon/Olympus/ Pilling/RFQ/Schöilly/Schott/ Storz/Stryker/Volpl/Winter&IBE/ Wolf/Zeiss	ACM/Aesculap/Dyonics/ Heine/HSW/Olympus/RFQ Storz/Stryker/Wolf/Zeiss

Function test

- Procedure:
- Hold one side of the fiberoptic light cable pointing toward a light source (e.g. window). For fiberoptic light cables, check the number of dark dots at the other end. These dark dots indicate broken optical fibers of the fibre optic bundle. Several broken optical fibres will not result in a noticeable impairment of the image quality. However, starting from about 20 to 25% of the light output surface, the light transmission and hence the image quality start to deteriorate considerably.
 - Check the fiberoptic light cable and the cable ends for damages such as burns or erosions.


- Information:
- Please note that the unscrewable adapters may potentially be re-used when ordering a new fiberoptic light cable.

1. Cleaning preparation

To achieve an optimal cleaning and disinfection outcome, both adapters should be disassembled from the fiberoptic light cable


Procedure Disassembly:

- To disassemble, screw the adapters from the respective end of the cable by turning anti clockwise.

Procedure Assembly:

- To assemble, screw the adapters onto the respective end of the cable by turning clockwise.

Information:

- Place the fiberoptic light cables onto a disposal screen in such a way that no rinsing-related shadows are created!
- Remove major contaminations from the fiberoptic light cable at the operating table.

General processing, information

The operator of the respective central sterile supply department (ZSVA / CSSD) as well as his/her assistants are responsible for the proper processing.

Medical devices intended for use in low-germ count or sterile conditions must be processed in compliance with the instructions from the manufacturer using adequate validated methods in such a way that the positive outcome of this method is clearly traceable and the safety and health of patients, users or third parties is not impaired.

Information:

- The manufacturer's specifications concerning concentration, temperature, reaction time etc. of the cleaning and disinfection agents used for the processing must be complied with.
- Only chemicals suitable and approved for processing optic fibre cables may be used.
- If the concentrations of the chemicals are too high, it may result in damages to the fibre optic cable and make the laser and etched lettering become illegible.
- The fibre optic cables should be processed promptly to prevent damage to optic fibre cables.
- Proper processing of the fiberoptic light cables starts at the operating table.
- New and repaired must pass through the entire processing cycle before being used!
- Fiberoptic light cables with visible damages must be repaired by an expert.
- Damaged fiberoptic light cables may impair the successful outcome of a surgical procedure.

2. Manual cleaning / disinfection

Auxiliary means:

- Synthetic brush with soft bristles

Procedure:

- Completely immerse the fiberoptic light cable into a cleaning and disinfection agent. Make sure that all fiberoptic light cables are fully immersed, i.e. that they are completely covered with cleaning and disinfection agent.
- Use a brush to clean the adapter threads on both ends of the fiberoptic light cables.
- Check the fiberoptic light cable for residual dirt and re-brush thoroughly, if necessary.
- After cleaning and disinfection, the fiberoptic light cable should be rinsed with fully desalinated water for the complete removal of any cleaning and disinfection agent residues.

Information:

- The respective instructions from the manufacturer concerning the concentration and reaction time of the cleaning and disinfection agent must be complied with.
- Adequate personal protection must be in place.
- Only brushes with soft synthetic bristles may be used to prevent damages to the fiberoptic light cable.
- The use of cleaning brushes made of metal and abrasive pellets (tipp cleaners) can lead to corrosion and damage the metal parts of the fiberoptic light cable.

3. Mechanical cleaning / disinfection

Use adequate validated process sequences for the cleaning and disinfection to demonstrably ensure a positive outcome (DIN EN ISO 15883-5). The individual process steps should be defined in the validation procedure by the respective central sterile supply department. To achieve an optimal cleaning result, the process parameters as well as the dosing quantity should be defined in co-operation with the manufacturer of the chemicals.

We recommend using the "Vario TD-programme" as cleaning and disinfection programme.

Auxiliary means:

- Washer-disinfector device

Procedure:

1. Place the fiberoptic light cables onto a disposal screen in such a way that no rinsing-related shadows are created.
2. Cleaning and disinfection
3. Visually check the fiberoptic light cables for cleanliness after they have been cleaned.


Fig. Schematic program run of the Vario TD program

Information:

- Fiberoptic light cables that are not completely clean must be re-cleaned (with manual pre-cleaning, if necessary).
 - To prevent damages to the fiberoptic light cable, it should be processed on a separate disposal screen in the cleaning and disinfection device.
 - If the pH for the cleaning is too high, it may cause damage to the fiberoptic light cables.
- The manufacturer of the Fiberoptic light cable recommends the following pH values for the cleaning:
- | | |
|----------------|--|
| pH 7.0 - pH 11 | Stainless steel (adapter & edging) |
| pH 4.0 - pH 11 | Plastics, silicone, glass surfaces (handle, tube & jacket) |

4. After-treatment / checking / packing

Auxiliary means:

- Magnifying glass
- Light source (eq. window)

Procedure:

- Before assembling the fiberoptic light cables, the individual parts must be checked for damages.
- The adapters should be re-mounted on the corresponding fiberoptic light cable before the sterilisation.
- Defective and non-functional fiberoptic light cables must be repaired or replaced.
- Use a standardised (DIN EN ISO 868/11607) approved packaging system to package the fiberoptic light cables.

Information:

- See also step 1, „Cleaning preparation“ and the point „Function test“

5. Sterilization

Adequate validated process sequences should be used for the sterilisation to demonstrably guarantee the positive outcome. We recommend using a fractionated vacuum procedure with saturated steam. Sterilisation with other sterilisation methods such as e.g. gas sterilisation or plasma sterilisation was not tested. If other sterilisation methods are used, the manufacturer of the fiberoptic light cable does not accept any liability in this context – the operator of the central sterile supply department will be solely liable.

Procedure, sterilization with saturated steam:

- Fractionated fore vacuum
- The sterilisation temperature must be between 134°C and 137°C with a stay-down time of 5 min. / 18 min. (vCJD) and a pressure of < 3050 mbar.

Information:

- Do not sterilise the fiberoptic light cables inside the supplied protective packaging.
- The operator of the respective central sterile supply department is responsible for the proper sterilisation.
- DIN EN ISO 17664 and DIN EN 285 must be complied with.
- National as well as international sterilisation guidelines must be complied with.
- The following temperatures must not be exceeded:
Stainless steel, Aluminium, Titanium alloys 150°C / Plastics 140°C

6. Disposal

No special precautions are required for the disposal of fiberoptic light cables.

General Information

Please always observe the applicable national regulations for the reprocessing of medical devices.

We will replace any demonstrably defective materials or deficient workmanship free of charge for the duration of two years from the invoice date. In warranty case, we are unable to cover the shipping costs or the risks associated with the shipping. Defective fiberoptic light cables may only be repaired by persons authorised by us and by means of original parts. Any repairs, alterations or upgrades not performed by the manufacturer or experts authorised by the manufacturer will void the warranty.

Disclaimer

The manufacturer of the fiberoptic light cable does not accept any liability for direct or consequential damages caused by the improper use, handling or improper processing, sterilisation and maintenance. Failure to follow the instructions, improper handling or use of products supplied by the manufacturer of the fiberoptic light cable for unintended purposes will completely void the warranty. The manufacturer of the fiberoptic light cable does not accept any liability whatsoever for damages arising hereof.

Status of Information

These operating instructions were issued in accordance with up-to-date facts and the date of revision is: **November 2014**

Manufacturer / Distribution

RfQ-Medizintechnik GmbH & Co. KG
Bruderhofstrasse 10-12
D-78532 Tuttlingen
Germany
Tel: +49 (0) 7461/9617-0
Fax: +49 (0) 7461/9617-20
e-mail: info@rfq.de