- DE Gebrauchsanweisung**
Starre medizinische Endoskope
- EN Instructions for use**
Rigid medical endoscopes
- FR Mode d'emploi**
Endoscopes médicaux rigides
- IT Istruzioni per l'uso**
Endoscopi rigidi ad uso medico
- ES Instrucciones para el uso**
Endoscopios rígidos médicos
- PT Instruções de utilização**
Endoscópios médicos rígidos
- NL Gebruiksaanwijzing**
Starre medische endoscoop
- PL Instrukcja obsługi**
Endoskopy medyczny sztywne
- ZH I 使用说明书**
医用硬管内窥镜

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Table of contents

1.	About this document	14
1.1	Purpose	14
1.2	Symbols used	14
2.	Intended use	15
3.	Safety information	15
4.	Testing, handling and maintenance	16
5.	Guideline conformity	16
6.	Description	17
6.1	Construction	17
6.2	Markings	17
6.3	Available designs and sizes	17
6.4	Combinable products	17
7.	Preparation for use	18
7.1	Visual inspection and function check	18
7.2	Provisioning	18
8.	Use	19
9.	Processing	19
9.1	Safe storage and transport	19
9.2	Cleaning and disinfection	19
	Manual cleaning / pre-cleaning and chemical disinfection	19
	Machine cleaning and thermal disinfection	20
	Removing deposits from optical end surfaces	21
9.3	Sterilisation	21
	Steam sterilisation (autoclaving)	21
	Fractionated pre-vacuum method	21
	Gravitation method	22
	Hydrogen peroxide sterilisation (STERRAD® method)	22
	Ethylene oxide sterilisation	22
9.4	Special precautions: Pathogens of Transmissible Spongiform Encephalopathy	22
9.5	Processing restrictions	22
10.	Assembly	23
11.	Disassembly	24
12.	Storage	24
13.	Service and maintenance	24
14.	Accessories / spare parts	25
15.	Disposal	25

1. About this document

1.1 Purpose

This document describes the correct handling and function of the rigid endoscope, as well as recommended processing methods.

This document may not be used to carry out endoscopic examinations or surgeries, nor may it be used for training purposes.

The respective current version of this document can be requested from Henke-Sass, Wolf.

If you as the user of this endoscope believe that you require more detailed information regarding the product's use and maintenance, please contact your representative.

1.2 Symbols used

The following symbols are used to make it easier for you to access the information:

	Instructions for preventing personal injury
	Instructions for preventing material damage
	Information to facilitate understanding or workflow optimisation
	Prerequisite
	Instruction

2. Intended use

Rigid medical endoscopes are used to visualise body cavities. Each endoscope was developed for diagnostic and surgical procedures in one of the following fields of application:

- Arthroscope: arthroscopic procedures
- Sinuscope: sinusoscopic procedures
- Otoscope: otoscopic procedures
- Laryngoscope: laryngoscopic procedures
- Laparoscope: laparoscopic procedures
- Hysteroscope: hysteroscopic procedures
- Cystoscope: cystoscopic procedures

For the benefit and safety of patients, physicians must select a method which they consider suitable based on their experience. If you as the user of this endoscope believe that you require more detailed information regarding the product's use and maintenance, please contact your representative.

3. Safety information

The endoscope may only be used by trained medical professionals in medical facilities.

- After delivery, inspect the endoscope for completeness and damage.
- Read, observe and keep the instructions for use.
- Use the endoscope only as intended, see "Intended use" on page 15..

For storage, transport and processing, ensure that the endoscope is not subjected to mechanical strain, particularly to prevent damage to the sensitive lens system.

WARNING

Risk of infection to the patient or medical professionals!

The endoscopes are delivered non-sterile as reusable products.

The state of the art and national laws require the observance of validated processes.

In general, users are responsible for the validation of their processes.

- **Ensure that the processing, material and personnel are suitable for achieving the results necessary.**
- **Observe any valid local operator regulations for all manual cleaning and drying processes.**
- **Clean / disinfect and sterilise the endoscope prior to initial use as well as each subsequent use of the endoscope.**
- **Bring the endoscope to the decontamination area after use. Observe valid protective measures to prevent contaminating the environment.**

WARNING

Risk of burns!

The optical fibres emit high-energy light at the distal end of the endoscope. This can cause the temperature of the body tissue to rise to 41 °C.

- **Avoid direct contact of the distal end with body tissue or flammable materials as it can cause burns.**
- **Reduce the light intensity of the cold light source when working near body tissue or flammable materials.**

WARNING

Risk of injury due to faulty endoscopes!

- **Carry out visual inspection and function check prior to each use.**
- **Only use endoscopes which are in perfect condition.**

4. Testing, handling and maintenance

Endoscopes from Henke-Sass, Wolf are precision medical instruments, and handling them requires great care.


- Inspect the endoscope for damage prior to and after use.
- If the endoscope is damaged, discontinue use and contact the manufacturer.
- Do not subject the endoscope to impact.
- Put the endoscope down carefully.
- Hold endoscope only by the ocular funnel / main part and not by the sheath.
- Do not bend the sheath.
- Do not bend the body after inserting the endoscope into the body. A piece broken off the endoscope can become lodged in the soft tissue or no longer appear in the endoscope's field of vision and thus remain in the body.
- Transport endoscopes individually and store them safely by using a screen basket or container.

5. Guideline conformity

<div>CE</div> <div>CE 0086</div>	The CE marking of the medical product complies with the guideline 93/42/EWG. It applies only when the product and/or packaging features this marking. In addition, products of Class IIa and higher are labelled with the identification number of the respective notified body.
For some repairs	
<div>CE 0535</div> <div>CE 0123</div>	

6. Description

6.1 Construction

- 1** Distal end
- 2** Sheath
- 3** Main part
- 4** Ocular funnel
- 5** Proximal end
- 6** Irradiation surface of the illumination fibres
- 7** Connection for illumination fibre, type ACMI
- 8** Adapter for illumination fibre, type Wolf, pre-assembled
- 9** Adapter for illumination fibre, type Storz / Olympus (assembly, see "Assembly" on page 23.)
- 10** For endoscopes with working channel: Adapter with stopcock (disassembly, see "Disassembly" on page 24.)

6.2 Markings

- Article number
- Serial number
- CE mark with identification number of the notified body where applicable: Endoscope conforms to the requirements of the guideline 93/42/EWG
- For autoclavable endoscopes: Writing "autoclavable"
- Specification of the direction of view
- Writing GERMANY
- Writing HENKE-SASS, WOLF
- For large-screen lenses: Writing **HM**
- For high-resolution optics developed specifically for the latest Full HD camera generations:
Writing **HD**
- For endoscopes of the third Henke-Sass, Wolf HD generation:
Writing **HD³**

6.3 Available designs and sizes

The endoscopes are available in the following designs and sizes:

- Straight endoscopes
- Angled endoscopes
- Endoscopes with working channel
- Sheath diameter 1.9–11 mm
- Working channel diameter 5.1–6.4 mm

6.4 Combinable products

You can combine the endoscopes with existing camera systems and with illumination fibres and instruments from Henke-Sass, Wolf.

7. Preparation for use

7.1 Visual inspection and function check

WARNING

Risk of injury due to faulty endoscopes!

- Carry out a visual inspection and function check prior to initial use as well as each additional use.
- Only use endoscopes which are in perfect condition.
- Clean / disinfect and sterilise the endoscope prior to initial use as well as each additional use of the endoscope. Contaminants on the irradiation surface of the illumination fibres can burn in during use, which impacts image quality.
- Ensure that the proximal end of the endoscope is dry to prevent the endoscope from fogging up during the examination / procedure.
- Ensure that no parts are missing or loose.
- Ensure that there are no residual cleaning agents or disinfectants on the endoscope.
- Inspect the entire endoscope, particularly the sheath, for contaminants and damage of any type, such as dents, scratches, cracks, bending and sharp edges.
- Inspect distal end, proximal end and irradiation surface of the illumination fibres for contamination and scratches. Make contaminants and scratches visible using light reflexes. Hold the connection of the optical fibres against the light and inspect whether the optical fibres illuminate evenly at the distal end.
- Check image quality: The image may not be blurry, clouded or dark. If necessary, remove deposits on the optical end surface using polishing paste provided, see “Removing deposits from optical end surfaces” on page 21..
- For endoscopes with locking device: Inspect between the sheath and the main part for contaminants and damage to ensure a fixed and secure connection.
- For endoscopes with working channel: Inspect all components of the adapter with stopcock for function and damage.
- Inspect free passage of the working channels.

7.2 Provisioning

- Clean / disinfect and sterilise the endoscope prior to initial use as well as each additional use of the endoscope, see “Processing” on page 19..
- Ensure that the proximal end of the endoscope is dry to prevent the endoscope from fogging up during the examination / procedure.
- If necessary mount adapter for illumination fibre, see “Assembly” on page 23..
- Mount illumination fibre (see manufacturer's specifications).
- If required, adapt the camera (see manufacturer's specifications).

8. Use

WARNING

Risk of burns!

The optical fibres emit high-energy light at the distal end of the endoscope. This can cause the temperature of the body tissue to rise to 41 °C.

- Avoid direct contact of the distal end with body tissue or flammable materials as it can cause burns.
- Reduce the light intensity of the cold light source when working near body tissue or flammable materials.

For endoscopes with working channel:

- If an instrument is inserted, do not move the lever on the stopcock because otherwise the instrument could be damaged or sheared off.
- Prepare the endoscope for processing immediately after use to prevent surface drying of the contaminants.

9. Processing

9.1 Safe storage and transport

If possible, reprocessing endoscopes immediately after use is recommended. Endoscope containers and trays are reusable. Trays must be inspected for visible contamination and cleaned prior to use. They can be cleaned manually or in an automatic cleaning unit using a cleaning agent.

- Always store endoscope securely and transport it to processing in a closed container to prevent damage to the endoscope and contamination of the environment.

9.2 Cleaning and disinfection

Manual cleaning / pre-cleaning and chemical disinfection

- Do not use fixating cleaning agents or hot water (>40 °C) as it can cause fixation of the contaminants and jeopardise successful cleaning.
- Do not scratch contaminants off with hard objects as this can cause damage to the optical end surfaces.
- Do not clean endoscope in an ultrasonic bath.

- ✓ Existing adapters are dismantled.
- ✓ For endoscopes with working channel: Adapter with stopcock is dismantled.
- ✓ For laryngoscopes: Handle is dismantled.
- Remove coarse contamination from the endoscope. Clean the endoscope with a soft brush under cold tap water until all visible contaminants have been removed.
- For working channels, holes and threads: Rinse with a water pistol for at least 10 seconds at a pressure of at least 3.8 bar (absolute).

The endoscopes have material compatibility with the Steris™ System 1 process.

- Disinfect endoscope. In doing so, observe the specifications of the disinfectant solution manufacturer regarding temperature, concentration and application time.

Non-compliance with the manufacturer's specifications can result in damage to the endoscope.

- Rinse endoscope with running water.
- If necessary, rinse working channel with a water pistol for at least 10 seconds at a pressure of at least 3.8 bar (absolute).
- If necessary, dry working channel with compressed air.
- Dry endoscope with a soft cloth.
- Carry out visual inspection, function check and servicing, see "Testing, handling and maintenance" on page 16..

Machine cleaning and thermal disinfection

The rigid endoscopes from Henke-Sass, Wolf are suitable for prevalent machine methods of cleaning and thermal disinfection. In doing so, use gentle cycles for rigid endoscopes and suitable cleaning agents and disinfectants. The instructions of the machine, cleaning agent and disinfectant manufacturers must be observed. The cleaning and disinfectant result must be confirmed by the machine, cleaning agent and disinfectant manufacturers in cooperation with the user.

The following methods have been validated by Henke-Sass, Wolf for the following rigid endoscopes:

- without working channel
 - with working channel (dia. ≥ 5 mm)
- Fix the endoscope to the loading rack in such a way that damage is prevented during cleaning.

The following materials and machines were used for the validation:


- Cleaning agent:
 - Alkaline: Neodisher FA; Dr. Weigert; Hamburg
 - Enzymatic: Endozime, Ruhof
- Neutraliser:
 - Neodisher Z; Dr. Weigert, Hamburg
- Cleaning and disinfecting unit:
 - Miele G 7736 CD
- Loading rack:
 - Loading rack E 327-06
 - MIC rack E 450
- Start cleaning process:
 - Pre-rinse with cold water for 1 minute
 - Drain
 - Pre-rinse with cold water for 3 minutes
 - Drain
 - Clean with 0.5% alkaline cleaning agent for 5 minutes at 55 °C or with 0.5% enzymatic cleaning agent at 45 °C
 - Drain
 - Neutralise for 3 minutes with warm tap water (<40 °C) and neutraliser
 - Drain
 - Intermediate rinse for 2 minutes with warm tap water (<40 °C)
 - Drain
- Carry out machine thermal disinfection considering the national requirements regarding the A_0 value (see DIN EN ISO 15883).
- Ensure that the exteriors of the endoscope are dry. If necessary, dry with a soft cloth.
- If necessary, dry working channels with compressed air.
- Carry out visual inspection, function check and servicing, see “Testing, handling and maintenance” on page 16..

Removing deposits from optical end surfaces

If deposits are found when checking the image quality, they can be removed with the provided polishing paste as follows:

i ➤ **Only clean with polishing paste if the image which you see through the endoscope is cloudy and blurry.**

- Apply polishing paste to a clean cotton swab.
- For large end surfaces: press cotton swab lightly on the end surface to be cleaned and rub it over the glass.
- For small end surfaces: place cotton swab lightly on the end surface to be cleaned and turn it.

- Clean all optical end surfaces with warm water and mild detergent to remove all polishing paste residue.
- Rinse optical end surfaces under running water.
- Dry optical end surfaces with a soft cloth.
- Clean / disinfect and sterilise the endoscope.
- Carry out visual inspection. If the deposits were not removed: send in endoscope for repair.

9.3 Sterilisation

i **Prior to each sterilisation, rigid endoscopes must be cleaned and disinfected according to the methods in these instructions for use.**

- Sterilise endoscopes in suitable packaging to prevent subsequent contamination.

Steam sterilisation (autoclaving)

In general, users are responsible for the validation of their processes.

Only endoscopes which are marked with the writing "autoclavable" are intended for autoclaving. The permissible processing methods are explained in the instructions at hand.

- When selecting the processing method, observe the valid national hygienic regulations and local provisions for hospital hygiene.

➤ **Comply with specified process parameters. The specified parameters have been validated to ensure the sterility of the endoscopes. Deviating process parameters could damage the endoscope. In this case, the guarantee and warranty shall become void.**

Autoclavable endoscopes can be sterilised with the French cycle (134 °C, 18 minutes, 3.1 bar (absolute) without restrictions regarding material compatibility.

- ✓ Existing adapters are dismantled.
- ✓ For endoscopes with working channel: Adapter for stopcock is dismantled.
- Sterilise endoscopes.
- When the sterilisation process has ended, allow the endoscopes to cool gradually to room temperature.

Fractionated pre-vacuum method

The following process has been validated:

Temperature	132–137 °C (270–278 °F)
Time	At least 3 minutes
Configuration	Double packed in sterilisation bags
Drying	At least 10 minutes

Gravitation method

The endoscopes have material compatibility with the gravitation method for a hold time of 15 minutes.

Hydrogen peroxide sterilisation (STERRAD® method)

Endoscopes without working channel and with working channel with a dia. ≥ 5 mm can be sterilised with the following STERRAD systems:

- STERRAD 100S
- STERRAD NX
- STERRAD 100NX

➤ Observe specifications of the manufacturer (ASP) regarding the corresponding method.

Ethylene oxide sterilisation

The endoscopes are material compatible with ethylene oxide sterilisation.

9.4 Special precautions: Pathogens of Transmissible Spongiform Encephalopathy

A comprehensive explanation of the necessary preventative measures with regard to agents of Transmissible Spongiform Encephalopathy (TSE) would go beyond the scope of this document.

It is assumed that pathogens of the Creutzfeldt Jakob Disease cannot be killed using normal disinfection and sterilisation processes. Therefore, the standard methods for decontamination and sterilisation are not sufficient if there is a risk of transferring Creutzfeldt Jakob Disease.

In general, only tissue with a low potential of TSE infection comes into contact with surgical instruments. In spite of this, special preventative measures must be taken for instruments which are used to treat patients with a known or suspected infection of TSE, as well as for patients at risk.

9.5 Processing restrictions


Repeated processing has only minimal effect on the endoscopes. The service life of the units is usually determined by wear and damage.

The endoscope can be damaged if the manufacturer's specifications are not observed.

➤ **Do not clean endoscope in an ultrasonic bath.**


10. Assembly

- 1** Connection for illumination fibre, type ACMI
- 2** Adapter type Wolf
- 3** Adapter type Storz / Olympus

- If necessary mount adapter for illumination fibre
- Ensure that the irradiation surface of the illumination fibre is clean.
- Mount illumination fibre (see manufacturer's specifications).
- If required, adapt the camera (see manufacturer's specifications).

For endoscopes with working channel:

- a** Union nut
- b** Stopcock
- c** Sealing cap, green
- d** Stopcock nut

- **In order to ensure sterility, only use grease which is suitable for medical instruments for the stopcock.**

- Lubricate stopcock **b**.
- Mount stopcock **b** and fix with stopcock nut **d**.
- Remove excess grease.

11. Disassembly

WARNING

Risk of burns!

Allow the illumination fibre to cool sufficiently before removing it. The ends can become very hot and cause serious burns.

➤ **Remove illumination fibre.**

➤ **Do not remove the ocular funnel because otherwise the endoscope will be damaged.**

- Unscrew existing adapters.
- For endoscopes with working channel, disassemble the adapter with stopcock.
 - Remove sealing cap **c**.
 - Loosen union nut **a**.
 - Loosen stopcock nut **d**.
 - Disassemble stopcock **b**.

12. Storage

Unsterile metal units must be stored in a clean, dry environment. The storage time of unsterile units is not limited; the units are made of a non-degradable material which maintains its stability when stored under the recommended conditions.

As long as endoscopes are stored unsterile in the original packaging, the following storage conditions apply:

- Temperature: $-10\text{ }^{\circ}\text{C}$ to $+40\text{ }^{\circ}\text{C}$
- Humidity: 10% to 90%
- Avoid direct sunlight.
- Store endoscope either in the original packaging or in a screen tray/container.
- Ensure that the endoscope is stored securely.
- Observe the respective valid national provisions when storing in a sterile condition.

13. Service and maintenance

Henke-Sass, Wolf does not supply original parts to independent workshops or other endoscope manufacturers. Thus only Henke-Sass, Wolf is in a position to carry out repairs using original parts. The original technical specifications and the operational safety of the endoscope can only be guaranteed by using original parts. The warranty for Henke-Sass, Wolf products shall become void if repairs are carried out by an workshop not authorised by Henke-Sass, Wolf. In this case Henke-Sass, Wolf is also no longer responsible for the technical specifications or safety of the product.

- Have the endoscopes repaired by Henke-Sass, Wolf only. For service, send the defective endoscope to the address of the sales partner.
- Clean, disinfect and sterilise the endoscope thoroughly prior to returning it for repair.
- Ideally, send in the endoscope in its original packaging. If this is not possible, package the endoscope to secure it for transport.

Henke-Sass, Wolf is not liable for damage resulting from improper shipping.

14. Accessories / spare parts

Designation	Article number
Polishing paste	8101001000
Adapter type Wolf	89530738
Adapter type Storz / Olympus	89542801
Sealing cap, green for 5 mm instruments (4 mm outlet)	89539010
Stopcock for endoscopes with working channel (article no. 7198004500 and 7198009000)	89509639
Stopcock for endoscopes with working channel (article no. 7198004511 and 7198009011)	89536169
Stopcock nut	89533648

15. Disposal

- Observe country-specific regulations and laws for the disposal of medical products.

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© copyright Henke-Sass, Wolf
Henke-Sass, Wolf GmbH