delmont
imaging

EN - Instructions for use

Universal light guide

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Guide de lumière universel

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Cable de fibra óptica universal

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TABLE OF CONTENTS

1. Intended use.....	3
2. Important information.....	3
3. Product description.....	3
4. Instructions on use and reprocessing.....	4
4.1. General principles	
4.2. Cleaning and disinfection	
4.3. Inspection and maintenance	
4.4. Packaging	
4.5. Sterilisation	
4.6. Storage	
4.7. Material stability	
4.8. Reusability	
4.9. Information on disposal	
5. Safety information.....	8

REF

1. Intended use

The universal light guide is used for endoscopic applications to transmit light during endoscopic diagnosis and therapy, or for other medical illumination applications when the properties and operating instructions are complied with.

The light guides are designed for use with halogen, xenon or LED based sources of cold light that are utilised in medical applications that involve endoscopes, medical instruments or microscopes and are equipped with an optical injection element for light guides with a bundle diameter of 4,8mm to 5,0mm or 3,5mm. High intensity visible light is conducted to the instrument through flexible light guides.

The use of the light guide together with LASER light sources and/or HF devices is not permitted. When selecting the light sources, follow the information provided under item 4, «Instructions on use and reprocessing», regarding the temperatures of use.

Contraindication: Contraindications related directly to the product are currently not known.

UMDNS code: 15-308: «Cable, fibre optic».

2. Important information

Symbols used:

Symbol	Meaning
	Caution
	Read the operating instructions prior to use
	Storage conditions (-20°C to +60°C)

3. Product description

The light guide is composed of an environmentally friendly SCHOTT PURAVIS® fibre optic bundle (type GOF70), which is hot fused to the light input surface. The end sleeves make it possible to connect any ordinary adapter sleeves (not included) which can be screwed on by the medical staff for all common endoscopes and light sources.

The light guide is distinguished by the following optical characteristics:

- Angle of acceptance (2α): $\geq 70^\circ$ at 587nm ($V_{(\lambda)}$, 1m length)
-> $NA \geq 0,57$
- Transmission : $\geq 60\%$ at 546nm (typically $\geq 65\%$)

4. Instructions on use and reprocessing

4.1. General principles

The light guide must be cleaned, disinfected and sterilised prior to each application. This is of particular importance when the light guide is used for the first time after delivery as all light guides are shipped in a non-sterile state (clean and disinfect after removing the transport packaging; sterilise after packaging). Efficient cleaning and disinfection is an essential requirement for effective sterilisation.

The user's responsibility with regard to ensuring sterility comprises:

- the exclusive use of device and product specific procedures that have been validated as sufficient for cleaning / disinfecting and sterilising,
- the assurance that all devices used (disinfector, steriliser) are serviced and checked on a regular basis,
- the assurance that the validated parameters are observed during every cycle.

The user is also required to observe all applicable statutory regulations and the rules established by the doctor's office or the hospital that concern hygiene. This applies in particular to the different regulations regarding an effective inactivation of prions.

4.2. Cleaning and disinfection

If possible, a mechanical procedure (disinfector) must be used for cleaning and disinfection. A manual procedure (even if involving the use of an ultrasonic bath) should not be used unless a mechanical procedure is not at hand due to the significantly lower effectiveness and reproducibility of a manual procedure.

The following pre-treatment must be performed in either case:

Coarse impurities must be removed from the light guide immediately following each application (within no more than 2 hours).

 **WARNING:** *Cleaning, disinfection and subsequent sterilisation with the adapter sleeves screwed on is not permitted. They must be removed and reprocessed separately prior to the reprocessing of the light guide.*

Use running water or a disinfectant solution for this purpose; the disinfectant should be aldehyde free (otherwise, blood contamination will become fixed in place), be of certified effectiveness (e.g. have VAH/DGHM or FDA approval or be labelled with the CE mark), be suitable for disinfecting the light guide and be compatible with the instruments (see chapter «Material stability»). Use only a soft brush and a soft cloth or cotton swabs to manually remove all visible impurities on the light guide and the light input and output surfaces, respectively. Use these items only for the purpose of removing impurities and refrain from using metal brushes or steel wool. Abrasive cloths or sponges may scratch these surfaces as well and are not suitable for cleaning.

Please note that the disinfectant used during the pre-treatment is only intended to provide operator protection and cannot substitute for the disinfecting step that needs to be completed later, when cleaning is finished.

Mechanical cleaning / disinfection:

When selecting the disinfector (combined cleaning and disinfecting unit), make sure:

- that the effectiveness of the disinfectant has been proven in principle (e.g. DGHM or FDA approval / CE mark in accordance with DIN EN ISO 15883),
- that, if possible, a certified program for thermal disinfection (A_0 value > 3000 or (for older units) at least 5min at 90°C) is employed (risk of disinfectant residue on the light guide during chemical disinfection),
- that the employed program is suitable for the light guide and comprised of a sufficient number of rinsing cycles,

- that you use only water that is sterile or has a low germ (max. 10 germs/mL) and a low endotoxin (max. 0,25 endotoxin units/mL) content (e.g. purified water / highly purified water),
- that the air used for drying is filtered,
- that the disinfectant is serviced and checked on a regular basis.

When selecting the cleaning agent system to be used, make sure:

- that this system is, in principle, suitable for cleaning instruments made of metals and plastics,
- that (if no thermal disinfection is employed) a suitable disinfectant with proven effectiveness (e.g. VAH/DGHM or FDA approval / CE mark) is used and that this disinfectant is compatible with the cleaning agent used,
- that the chemicals used are compatible with the light guides (see chapter «Material stability»).

The concentrations specified by the manufacturer of the cleaning agent and, if applicable, the disinfectant must be strictly observed.

The process is:

1. Place the light guide in the disinfectant. Make sure the instruments do not come into contact with one another and the minimum permissible bending radius is 50mm.
2. Start the program.
3. Remove the instruments from the disinfectant when the program is complete.
4. Check and package the light guides as quickly as possible after taking them out, if necessary after additional drying in a clean location.

The general suitability of the instruments for efficient mechanical cleaning and disinfection was proven by an independent, accredited test laboratory using the disinfectant G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent NEODISHER MediClean forte (Dr. Weigert GmbH & Co. KG, Hamburg). The procedure described above was applied for this process.

Manual cleaning and disinfection:

When selecting the cleaning agents and disinfectants to be used, make sure:

- that they are, in principle, suitable for cleaning / disinfecting instruments made of metals and plastics,
- that the cleaning agent (if applicable) is suitable for ultrasonic cleaning (no foam generation),
- that a suitable disinfectant with proven effectiveness (e.g. VAH/DGHM or FDA approval / CE mark) is used and that this disinfectant is compatible with the cleaning agent used,
- that the chemicals used are compatible with the light guides (see chapter «Material stability»).

Combined cleaning agents / disinfectants should be avoided if possible. Combined cleaning agents / disinfectants may only be used if contamination is low (no visible impurities).

The concentrations and residence times specified by the manufacturer of the cleaning agents and disinfectants as well as the specifications for rinsing must be strictly observed. Use only freshly made solutions as well as water that is sterile or has a low germ (max. 10 germs/mL) and a low endotoxin (max. 0,25 endotoxin units/mL) content (e.g. purified water / highly purified water) and filtered air for drying.

The process is:

Cleaning:

1. Immerse the light guide in the cleaning bath for the specified residence time. While doing so, make sure the light guides do not come into contact with one another and are covered sufficiently. Completely brush off the light guide, using a soft brush after half of the residence time has elapsed. **But do not brush the light input and output surfaces**, use only a soft cloth or cotton swabs for cleaning the light input and output surfaces.
2. Proceed by removing the light guides from the cleaning bath and rinse them completely at least three times for 1 minute under running water.

3. Check the light guides for any damage.

Disinfection:

4. Immerse the cleaned and checked light guide in the disinfectant bath for the specified residence time, making sure that the light guides are covered sufficiently. Ensure that the light guides are not in contact with one another.

5. Proceed by removing the light guides from the disinfectant bath and rinse them completely at least five times for 1 minute under running water.

6. Dry the light guides by blowing them off with filtered compressed air.

7. Package the instruments as quickly as possible after taking them out. If necessary after additional drying in a clean location.

The general suitability of the instruments for efficient mechanical cleaning and disinfection was proven by an independent, accredited test laboratory using the cleaning agent Cidezyme/Enzol and the disinfectant Cidex opa (Johnson & Johnson GmbH, Norderstedt). The procedure described above was applied for this process.

4.3. Inspection and maintenance

Inspection: When cleaning or cleaning / disinfection is complete, check all light guides for corrosion, damaged surfaces, chippings and impurities and sort out any damaged light guides. The light guides must be disposed of after no more than 100 reprocessing cycles (see chapter «Reusability»). If still soiled, the light guides must be cleaned and disinfected again.

Maintenance: Neither maintenance nor installation are required. The use of instrument oils is not permitted.

4.4. Packaging

Please package the light guide in single use packages (single or twin packaging) made of paper/film (mandatory for steam sterilisation) or Tyvek/film (mandatory for LTP sterilisation) which meet the following requirements:

- DIN EN ISO/ANSI AAMI ISO 11607.
- Suitable for steam sterilisation (min. temperature resistance up to 141°C (286°F), sufficient vapour permeability) or LTP sterilisation (Tyvek).
- Sufficient protection of the light guides or sterilisation packages against mechanical damage.

4.5. Sterilisation

Only the sterilisation procedures listed below may be used for sterilisation; other sterilisation procedures are not permissible:

Steam sterilisation:

- Fractionated vacuum process or gravitation procedure¹ (with sufficient product drying).
- Steam steriliser according to DIN EN 13060 or DIN EN 285.
- Validated in accordance with DIN EN ISO 17665 (formerly: DIN EN 554/ANSI AAMI ISO 11134) (valid IQ/OQ (commissioning) and product-specific performance qualification (PQ)).
- Maximum sterilisation temperature 138°C (280°F; plus tolerance according to DIN EN ISO 17665 (formerly: DIN EN 554/ANSI AAMI ISO 11134)).
- Sterilisation time (exposure time with sterilisation temperature) at least 5min² (gravitation procedure) and 3min² (fractionated vacuum process) at 132°C (270°F)/134°C (273°F).

¹ The use of the less effective gravitation procedure is not permitted unless the fractionated vacuum process is not available.

² Or 18 min (prion inactivation).

WARNING:

- As a basic principle, the use of Tyvek/film sterilisation packaging is not permissible for steam sterilisation; only the use of paper/film sterilisation packaging is permitted! The use of a sterilisation container requires that the client be responsible for conducting inspections that are specific to the product, load capacity, procedure, steriliser and establish that the sterilisation efficiency is unaffected.
- The use of flash sterilisation or other procedures that are intended exclusively for use without sterilisation packaging is not permissible under any circumstances!

The general suitability of the instruments for efficient steam sterilisation was proven by an independent, accredited test laboratory using the steam steriliser Systec V-150 (Systec GmbH Labor-Systemtechnik, Wettenberg); this proof has been established for both the fractionated vacuum process and the gravitation procedure. This process considered both conditions typically found in hospitals and doctor's offices and the procedure described above.

Low-temperature plasma (LTP) sterilisation:

- LTP steriliser in accordance with DIN EN ISO 14937.
- Validated according to DIN EN ISO 14937 (valid IQ/OQ (commissioning) and product-specific performance qualification (PQ)).
- Maximum sterilisation temperature 60°C (140°F; including permissible tolerance).
- Long cycle (Sterrad 100S).

WARNING:

- As a basic principle, the use of paper/film sterilisation packaging is not permissible for LTP sterilisation; only the use of Tyvek/film sterilisation packaging is permitted!
- The use of short cycles is not permissible under any circumstances!

The general suitability for efficient LTP sterilisation was proven by an independent, accredited test laboratory using the STERRAD 100S procedure (ASP Inc., U.S.A.) and the long cycle/a sterilisation period (H₂O₂ exposure time during sterilisation temperature plus plasma phase) of 20min at 45-55°C (113-131°F) and 6 mg/L H₂O₂. The procedure described above was applied for this process. When applying the LTP procedure, observe the applicable regulations in your country; certain application restrictions may apply.

Do not use heated air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation!

4.6. Storage

Once sterilised, the light guide must be stored dry, dust free and protected against light in the sterilisation packaging.

Storage temperature: -20°C to +60°C.

4.7. Material stability

When choosing the cleaning agent and disinfectant, make sure that they do not contain the following components:

- organic, mineral and oxidising acids (minimum permissible pH-value 5,5),
- strong bases (maximum permissible pH-value 11; neutral/enzymatic or weakly alkaline cleaning agents are recommended),
- organic solvents (such as alcohols, ether, ketone, and gasoline),
- oxidising agents (e.g. hydrogen peroxides),
- halogens (chlorine, iodine, and bromine),

- aromatic/halogenated hydrocarbons,
- tri-/perchloroethylene.

Never use metal brushes, steel wool, abrasive cloths or sponges to clean the light guides.

The light guides must never be exposed to temperatures exceeding 141°C (286°F)!

Temperatures of use:

For the overall light guide: +10°C to +40°C.

For the light input and light output surface: max. 350°C on the input side (long sleeve),
 max. 150°C on the output side (short sleeve).

4.8. Reusability

If handled with the necessary care and remaining undamaged and unsoiled, the light guides can be reused up to 100 times; using them beyond this limit or using damaged and/or soiled instruments is the responsibility of the user.

No liability in case of non-observance.

The fact that the light guide is reusable up to 100 times was validated in the following conditions:

Cleaning and disinfection: Mechanical cleaning and disinfection in MIELE Professional PG8535 disinfectant using NEODISHER MediClean forte (0,5 vol.%) cleaner; thermal disinfection at 95°C with a residence time of 10min.

Sterilisation: Fractionated vacuum process (triple) in a LAUTENSCHLÄGER PROTOCERT 839 steam steriliser (134°C, 10 min., 3 bar).

4.9. Information on disposal

After use and prior to disposal, the light guide must be disinfected or disposed of as hazardous medical waste.

5. Safety information

WARNING:

- Immediately before and after every use, the light guide and accessories used with it must be examined for fault free state, correct functioning, and completeness!
- Damaged light guides must be withdrawn from use immediately!
- Because of the high amounts of light energy, the light input surfaces and illuminated surfaces can heat up to an extent that may result in burns!
- Light sources may be used for the above purpose only if the maximum temperatures of use (see above) are complied with.
- Never look into the open end of a light guide which is connected to a light source! Risk of blinding!
- The light guide is not provided in a sterile state. Consequently, before first use and every subsequent use, the light guide must be cleaned, disinfected, and/or sterilised by means of validated processes!
- A minimum bending radius of 50mm is permissible in the flexible section of the hose.
- Installation information: the light guide may not be bent or crushed at any time.
- The light guide is not suitable for mechanical, dynamic movements. It may only be subjected to manual bending movements that are within the limits of the minimum bending radius.
- The light guide must not be used under tension during installation and use.
- The light guide is not intended for invasive use!

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