delmont
imaging

EN - Instructions for use
iCare imaging station

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TABLE OF CONTENTS

1.	General presentation.....	4
2.	Safety instructions.....	4
2.1.	General instructions	
2.2.	Instructions specific to the light source	
2.3.	Contraindication	
3.	Regulatory notice.....	5
3.1.	Compliance	
3.2.	Electromagnetic interferences and electrostatic discharges	
3.3.	Materiovigilance	
3.4.	Life-end	
4.	Product description.....	6
4.1.	Control unit description	
4.2.	Camera head description	
5.	Installation.....	8
6.	Functioning.....	9
6.1.	Powering on	
6.2.	Connecting the endoscope to the camera head	
6.3.	White balance	
6.4.	Focus	
6.5.	Other camera head buttons features	
6.6.	Turning off the product	
7.	Disinfection procedure of the camera head.....	10
8.	Failure detection.....	10
8.1.	The light indicator of the front panel does not illuminate when power is turned on	
8.2.	The LED source lights up but the luminous flux is insufficient	
8.3.	There is no more light but the fans continue to work	
8.4.	The light indicator of the front panel is fixed but no image appears on the screen	
8.5.	The image is blurred	
9.	After-sales service and maintenance.....	11
10.	Technical data.....	11

TABLE OF CONTENTS

11. Electromagnetic compatibility.....	12
11.1. Guidance and manufacturer's declaration: electromagnetic emissions	
11.2. Guidance and manufacturer's declaration: electromagnetic immunity	
11.3. Recommended distances between portable and mobile RF communication systems for this product	
12. Used symbols.....	15

1. General presentation

iCare is a combined HD CMOS camera and LED light source in a single control unit, designed for use by qualified gynecologists when performing surgical endoscopy or diagnostic procedures.

It includes all the following elements:

- a control unit that contains a 64 Watts LED light source,
- an HD CMOS sensor with its 22mm C-mount focal length lens, 3 pre-programmed buttons, 2.99mm cable, connector and sealing cap,
- an EU power cable,
- an HDMI cable,
- a USB key,
- two WiFi keys that allow our imagyn application to work properly (for more information on the imagyn software, please refer to the corresponding instructions for use).

On the camera side, it is a 1/4" colour CMOS mono HD camera with remote electronics. Its small ergonomic sensor, automatic shutter, sensitivity, excellent resolution and exceptional colour rendering make it the ideal medical tool.

On the light source side, this LED source is specially designed for use in various endoscopic diagnostic or surgical applications. Its ease of use and powerful lighting make it the ideal multidisciplinary medical tool.

To get the most out of iCare while taking all the necessary precautions, it is essential that you read this manual.

Phrases with the  symbol correspond to points that require special attention.

Phrases with the  symbol are information.

 *This equipment was delivered to you in a cardboard box. Keep it for possible transport.*

2. Safety instructions

2.1. General instructions

- Read the instructions for use.
 - Comply with the conditions of use and storage.
 - The unit may only be opened by a competent technician authorised by the manufacturer.
 - Do not insert metal objects into the unit to avoid electric shock, fire, short circuit or dangerous emission.
 - Do not expose the unit to water splashes or in a place that is too humid.
 - Use only accessories supplied with the unit or offered as an option by the manufacturer.
 - Do not place heavy objects on the unit.
 - This unit is not sterile.
 - If the power cord is damaged, turn off the power immediately. It is dangerous to operate this unit with a damaged cord. To disconnect the cord, pull it by the plug. Never pull on the cord itself.
 - Unplug the unit from power if it is not going to be used for a few days or more.
 - Ensure sufficient air circulation to avoid overheating inside the unit: at least 15cm all around the unit.
- Do not cover it and make sure that the feet of the unit are present.

- Do not use corrosive or abrasive products to clean the device but only the disinfectant liquids recommended in the corresponding chapter.
- Before each use, make sure that the unit has no rough surfaces, sharp edges or protrusions that could cause safety problems.
- The surface temperature can exceed 41°C (after a few minutes of use). Therefore, avoid contact with skin.
- To avoid the risk of electric shock, this device may only be connected to a power supply equipped with a protective ground.
- The devices connected to the inputs/outputs must comply with IEC 60950-1.
- This device may not be modified without the manufacturer's permission. If the medical device is modified, a check and test must be carried out to ensure that the medical device complies with the safety instructions.
- Check with the manufacturer the compatibility of your endoscope before use.
- To prevent errors or delays in diagnosis, it is important to ensure that the monitor settings used are optimized for the procedure being performed so that a clear, noise-free colour image is obtained.
- This device should be used on individuals (patients) suitable for endoscopic procedures.
- Applied parts of electro-medical devices that can be used in conjunction with iCare shall be of BF-type or CF according to standard 60601-2-18. Check this compatibility before each operation for safe use.

 *It is recommended to have a second camera and light source available to use in case an absence or degradation of performance is observed.*

2.2. Instructions specific to the light source

- Never look at the light output or the end of the light cable.
- Do not insert anything other than a light cable into the hole provided for this purpose, otherwise the optical system may be damaged.
- Use of cables or accessories other than those specified may result in increased emissions or decreased immunity of the unit.
- After using the source, when removing the fiber from its light guide, the temperature on the metal coupler of the fiber is very high and can cause burns.
- The light power at the light cable outlet can cause eye damage. Handle the light cable carefully when the device is in use.
- Do not place the distal end of the light cable or endoscope directly on the patient or on any other flammable material (sheets, gauze, operating fields, etc.) as this can be very hot and cause burns.

 *This product is equipped with Group 1 LEDs according to IEC 62471. Do not look directly at the light to avoid any ocular risk.*

2.3. Contraindication

The use of iCare is contraindicated when endoscopic practice is contraindicated for the patient.

3. Regulatory notice

3.1. Compliance

This product was designed and manufactured by a company with a certified quality system.

It meets the requirements of European Directive 93/42/EEC on medical devices. It therefore complies with the relevant electrical safety (IEC) and electromagnetic compatibility (EMC) standards.

3.2. Electromagnetic interferences and electrostatic discharges

Although this product complies with EMC standards, it is possible that under very special circumstances it may cause interference to other devices, or may itself be affected by other devices or an adverse electromagnetic environment.

To avoid these situations, it is recommended:

- to ensure the quality of the electrical network (especially the grounding of all equipment and trolleys),
- to keep the device away from electromagnetic sources (for example, a compressor, a motor, a transformer, an HF generator, etc.).

3.3. Materiovigilance

As any medical device, this device is subject to the provisions of the materiovigilance. Any serious malfunction must therefore be reported to the competent authorities and the manufacturer as soon as possible and with the greatest possible precision.

Manufacturer's contact information: refer to the last page of this IFU.

3.4. Life-end

This device carries the recycling symbol in accordance with the European Directive 2002/96/EC on the Waste of Electrical and Electronic Equipment (WEEE).

Proper disposal of this device will help prevent any harmful consequences for the environment and human health. The  symbol on the unit or accompanying documentation indicates that this product may under no circumstances be treated as household waste. It must therefore be handed over to a waste collection centre responsible for recycling electrical and electronic equipment.

For disposal, observe the waste disposal regulations in force in the country of installation.

For further details on the processing, recovery and recycling of this unit, please contact your nearest commercial representative for instructions.

4. Product description

4.1. Control unit description

[1]: The fans on the rear of the unit must not be blocked to prevent overheating. The light source is equipped with an automatic safety feature that stops the lighting if the internal temperature becomes excessive. Fans on the left side and under the unit are also present for better ventilation.

[2] : Power is supplied through the rear panel power outlet, which must be connected to the power supply via the cord provided with the product. This socket has a fuse cover and the main switch for powering on. When replacing a fuse, it is imperative to disconnect the product from the network and use a fuse of the same type. The T of «T2A» means «time-delay». Use only fuses marked UL/CSA.

[3] : The unit has an HDMI video output.

Devices that connect to the «VIDEO OUT» sockets must comply with IEC 60950.

[4] : The unit has a USB output for a storage key.

[5] : The unit has two USB outputs for two WiFi keys.

[6] : The unit has an equipotential plug, which can be connected to other electromedical devices in order to reduce the formation of different electrical potentials.

[7] : The indications and symbols on the rear of the unit enable it to be identified in accordance with international standards IEC 60601-1, IEC 60601-2-18, IEC 60417 and EN 980 (see the corresponding section 12).

[8] : The unit has a connector for plugging in the light cable.

[9] : The unit has a light indicator that alerts the user to the different phases of use of the product.

[10] : The unit has a connection for plugging in the iCare camera head.

[11] : The device has a stand-by button.

4.2. Camera head description

[12] : This button allows to launch a white balance by long press.

[13] : This button allows to switch from a minimum value of light power to a maximum value by short press, and vice versa. These values are 75% and 100% of the light output respectively.

[14] : This button allows:

- by short press, to launch an image capture,
- by long press, to start or stop a video recording.

5. Installation

No special training is required to install this medical device. Please refer to the instructions in this section of the user manual.

Once all iCare components have been unpacked from their original packaging, please perform the following actions:

- Place the control unit on a level and stable surface. If you place it in a compartment, make sure that it is sufficiently ventilated (at least 15cm around the unit).
- Connect the power cord to the power outlet [2] on the rear of the unit.
- Connect the other end of the power cord to an electrical outlet.
- Connect the HDMI cable to the corresponding output [3] on the rear of the unit.
- Connect the second end of the HDMI cable to the corresponding input on the monitor.
- Place the two WiFi keys on the two USB outputs [5] on the rear panel.
- Place the USB key on the USB output [4] on the rear panel.
- Set the power supply button on the rear panel of the unit to position « I ». The unit enters in a booting sequence: the indicator light [9] on the front panel flashes rapidly. This booting sequence ends after approximately 1 minute 30. Then, the unit goes into standby mode: the front panel indicator light flashes slowly.
- Power on the monitor. A colour bar should appear on the screen.
- Plug the camera head connector to the front of the control unit [10]. A red coding is present on the connector of the camera head and above the corresponding socket of the control unit. Align these two coded pins to connect the camera head. A « click » lock must be heard.
- Insert the light cable into the hole provided for this purpose [8]. A « click » lock must be heard. Connect the other end of the light cable to your endoscope.
- Place an endoscope in front of the camera lens.
- Press the stand-by button [11] on the right side of the control unit to start the device: the indicator light on the front panel becomes fixed. Another press will put it back on standby.

It is then possible to check the correct operation of the camera and especially to ensure that the settings of the monitor give full satisfaction thanks to the colour bar which appears when disconnecting the sensor.

6. Functioning

6.1. Powering on

The device is equipped with a switch placed on the rear of the device [2]. It is activated by changing this switch to position « I ». The light source is set to the last power used or to the value saved by the practitioner in imagyn.

6.2. Connecting the endoscope to the camera head

Bring the dewclaws on the ring of the lens closer together until the endoscope can be inserted. Once the endoscope is inserted, release the two dewclaws.

The endoscope can be unlocked by simply bringing the dewclaws closer together.

6.3. White balance

White balance is triggered by pressing and holding the [12] button on the camera head. The following procedure should be followed:

- Once the camera is paired with the endoscope and the light source is activated, film a white surface.
- Start the white balance by pressing the corresponding button.
- Continue filming the white surface as long as the message «Processing AWB...» is present on the screen.
- Keep the button pressed until the message «AWB OK» appears on the screen. This delay avoids unintentional press.

6.4. Focus

Use the lens ring associated to the camera head to focus.

Once the endoscope is connected and the light source activated, slowly turn the ring to find a position where the objects observed are sharp.

Focusing far enough away allows a sufficient depth of field to be obtained for the operation, thus avoiding regular focusing.

6.5. Other camera head buttons features

The [13] button on the camera head switches from a A value to a B value of the light source. These basic values are set at 75% and 100% light output but can be adjusted via imagyn.

The [14] button allows:

- by short press, to launch an image capture,
- by long press, to start a video recording. Another long press stops the video recording.

6.6. Turning off the product

To turn off the device, put the switch on the rear panel of the control unit [2] to the « 0 » position.

7. Disinfection procedure of the camera head

Pre-disinfection cleaning	Disinfection	Sterilization	
Immersion in an enzymatic cleaning solution (Aniosyme DD1, Hexanios G+R or surface cleaning (Anios wipe))	Immersion in a glutaraldehyde solution (Anios-Steranios Laboratory 2%)	Sterilization by STERRAD process	Sterilization with ethylene oxide

 Any other method of disinfection is prohibited. Damage caused by these other methods cannot be borne by the manufacturer.

 The use of a sterile protective cover on the camera is recommended regardless of the procedure followed in sterilization service and for the duration of the surgical procedure.

 The camera is not autoclavable.

 The camera is not compatible with automatic washers/disinfectors.

 Existing alkaline solutions for pre-disinfection of certain medical devices are FORBIDDEN for pre-disinfection of our cameras.

 Parts that have been in contact with the disinfectant must be rinsed thoroughly.

 Use non-woven compresses to dry the optics so as not to scratch them.

 The procedures described in this chapter are provided as advice, not in any way to replace official recommendations or guidelines.

 Decontamination being linked to the selected products, methods and/or tools, remains the sole responsibility of the personnel concerned.

8. Failure detection

8.1. The light indicator of the front panel does not illuminate when power is turned on

- Check that the power cord is connected to the network and the device, and that the power switch on the rear panel of the unit is in the « I » position.

- Check fuses for proper condition (use T2A - 250V - UL/CSA fuses only).

8.2. The LED source lights up but the luminous flux is insufficient

Check the actual value of the light output via imagyn. If necessary, put it at 100%. If the problem persists, check that the light cable is properly connected. If necessary, check the condition of your light cable and optics.

8.3. There is no more light but the fans continue to work

Check that the light cable is properly connected. If this is the case, wait a few minutes: the light source is equipped with a safety device that cuts off the power supply to the LED if the temperature inside the device is too high. Once the temperature has dropped, the source can be used again.

As a reminder, there must be enough space around the unit to ensure sufficient cooling (15cm all around) and nothing must obstruct the fan grids below, on the back and on the left side of the device.

If the defect persists and it is necessary to return the device to the after-sales service, take care to send it to us in its original packaging after having disinfected it.

8.4. The light indicator of the front panel is fixed but no image appears on the screen

- Check that the camera head is connected to the control unit (otherwise a colour bar will be displayed).
- Check that the control unit is correctly connected to the monitor (cable in good condition and plugs properly inserted).
- Check that the monitor is turned on, that the correct video input is selected and that the screen image settings are not in the minimum position (colour, light and contrast).
- Check for the presence of light by inspecting the light source, light cable and endoscope.

8.5. The image is blurred

- Check that there is no fog on the camera lens or endoscope.
- Check the focus of the lens.

9. After-sales service and maintenance

No specific maintenance is required on this device.

Simply make sure that the aerators are not clogged with dust. If this is the case, unplug the unit and vacuum the dust.

 *Errors of use are not covered by the warranty.*

If a defect persists and it is necessary to return the device to the after-sales service, take care to ship it in its original packaging. Also, the entire device (control unit, camera head, lens, cables, etc.) must be returned. Please attach a small explanatory note to the shipment regarding the defect found.

 *The equipment must be disinfected before returning for repair.*

When returning the equipment, check its condition and express reserves on the delivery order if necessary, confirming them to the carrier by registered letter within 48 hours. In the case where a material dispatched by us would suffer damage during its transport, the amount of the repair is charged either to the carrier if the reserves were made within the time limits, or to the recipient in the contrary case.

In case of an incident, please contact our customer service or our nearest sales representative.

10. Technical data

Camera side:

- Waterproof HD CMOS sensor (IP67)
- 22 mm C-mount lens
- Resolution: 1920 x 1080
- Interlaced scanning
- Definition > 900 lines
- Sensibility : 2 000 lux at F8
- Signal-to-noise ratio: 54 dB
- Automatic electronic shutter (1/50 to 1/10 000)

- Pre-programmed buttons on the camera head
- Colour bar
- White balance
- 1 HDMI output
- 2 USB outputs for WiFi connection
- 1 USB output for a storage key

Light source side:

- LED technology
- Nominal power: 64 W
- Colour temperature: 6 000°K
- LED typical lifetime: 50 000 hours
- Compatible light cable type: Storz
- Automatic thermal protection system
- Automatic detection of presence/absence of the light cable

General specifications:

- Power supply: 100-230 V ~; 50/60 Hz
- Power consumption: 130 VA
- Two fuses T2A - 250V UL/CSA 5 x 20 mm
- Equipotential plug
- Dimensions of the control unit (w x h x d): 310 x 75 x 310 mm
- Weight of the control unit: 4 200 g

Environment :

- Working temperature: +10°C/+40°C
- Working humidity: 30 to 75%
- Transport and storage temperature: -10°C/+40°C
- Transport and storage humidity: 20 to 85%
- Working, transport and storage atmospheric pressure: 800 hPa to 1 060 hPa
- Not protected against water splashes (IPX0)
- Not suitable for use in the presence of flammable anesthetic mixture with air with oxygen or nitrous oxide

Regulation :

- Electrical safety class 1, type BF
- Complies with European directive 93/42/EEC
- Complies with international norms: IEC 60601-1, IEC 60601-2-18, IEC 60601-1-2, IEC 60417 and EN 80
- IEC 62471: group 1 risk

11. Electromagnetic compatibility

11.1. Guidance and manufacturer's declaration: electromagnetic emissions

The « CMOS camera + LED source » reference equipment is intended for use in the electromagnetic environment specified below. The user must ensure that it is actually used in this environment.

Emissions test	Compliance	Electromagnetic environment: Guidance
RF emissions CISPR 11	Group 1	This « CMOS camera + LED source » product only uses radio power for its subsystems. It therefore emits very low RF energy and is not likely to interfere with nearby electronic devices. This « CMOS camera + LED source » product must be used in all installations, other than residential installations and premises directly connected to the public low voltage power distribution network intended to supply residential buildings.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Compliant	
Voltage fluctuations/Flicker IEC 61000-3-3	Compliant	

11.2. Guidance and manufacturer's declaration: electromagnetic immunity

The « CMOS camera + LED source » reference equipment has been designed for use in the electromagnetic environment specified below. The user must ensure that it is actually used in this environment.

Immunity test	IEC 60601 Severity Level	Compliance Level	Electromagnetic environment: Guidance
Electrostatic discharges IEC 61000-4-2	± 6 kV via contact ± 8 kV via air	± 6 kV ± 8 kV	The floor must be made of wood, concrete or tiles. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.
Rapid transient peaks IEC 61000-4-4	± 2 kV power lines ± 1 kV input/output lines	± 2 kV ± 1 kV	The quality of the main power supply must be the one of a typical commercial or hospital environment.
Electric shocks IEC 61000-4-5	Differential mode ± 1 kV Common mode ± 2 kV	± 1 kV ± 2 kV	
Power outages, short power outages and voltage fluctuations IEC 61000-4-11	<5% UT - for 10 ms 40% UT - for 100 ms 70% UT - for 500 ms <5% UT - for 5 ms	5% UT 10 ms <40% UT 100 ms <70% UT 500 ms <5% UT 5 s	The quality of the main power supply must be the one of a typical commercial or hospital environment. If the user of this product must be able to continue working during power outages, it is recommended that this product be powered by UPS or battery.
Magnetic field at mains frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The magnetic field at the mains frequency must be at a characteristic level of a location (50/60 Hz) in a typical commercial or hospital environment.

Note: UT is the nominal value of the electrical voltage applied during the test.

Guidance and manufacturer's declaration: electromagnetic emissions

The « CMOS camera + LED source » reference equipment is intended for use in the electromagnetic environment specified below. The user must ensure that it is actually used in this environment.

Safety test	IEC 60601 Severity Level	Compliance Level	Electromagnetic environment: Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	Portable and mobile RF communication devices should not be used at a distance, including cables, from this product that is less than the recommended distance, calculated by applying the formula that corresponds to the transmitter frequency. $d = 1,16 \sqrt{P}$ $d = 1,16 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,33 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ Where P is the maximum output power of the transmitter, in Watts (W), assigned by its manufacturer and (d) is the recommended separation distance in metres (m). Field strength levels emitted by fixed RF transmitters - which must be established by in situ electromagnetic measurement - must be below the compliance level in each frequency band. Interference may occur with devices marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the highest frequency band should be used.

Note 2: These recommendations may not be applicable in all situations. The propagation of electromagnetic waves is altered by absorption and reflection from structures, objects and people. For transmitters whose maximum output power is not listed in the table above, the recommended separation distance d, in metres (m) can be established using the equation applicable to the transmitter frequency, where P is the maximum output power of the transmitter in Watts (W) assigned by the transmitter manufacturer.

11.3. Recommended distances between portable and mobile RF communication systems for this product

The « CMOS camera + LED source » reference equipment has been designed for use in the electromagnetic environment in which the emitted RF interference is controlled. The user of this equipment can help avoid electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication systems (transmitters) and this equipment, as recommended below, as a function of the maximum output power of the communication system.

Maximum assigned transmitter output power in W	Separation distance as a function of transmitter frequency (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 1,16 \sqrt{P}$	$d = 1,16 \sqrt{P}$	$d = 2,33 \sqrt{P}$
0.01	0.116	0.116	0.233
0.1	0.366	0.366	0.736
1	1.16	1.16	2.33
10	3.66	3.66	7.36
100	11.6	11.6	23.3

Note 1: At 80 MHz and 800 MHz, the separation distance given in the upper frequency band applies.

Note 2: These recommendations may not be applicable in all situations. The propagation of electromagnetic waves is altered by absorption and reflection from structures, objects and people. For transmitters whose maximum output power is not listed in the table above, the recommended separation distance d , in metres (m) can be established using the equation applicable to the transmitter frequency, where P is the maximum output power of the transmitter in Watts (W) assigned by the transmitter manufacturer.

12. Used symbols

Symbol	Description		
	Symbol for «Manufacturer»	Legal Manufacturer: DELMONT IMAGING Zone Athélia V 390 Avenue du Mistral 13600 La Ciotat, FRANCE	Tel.: +33 (0) 9 51 51 30 30 Fax: +33 (0) 9 57 51 31 00 Email: contact@delmont-imaging.com Website: www.delmont-imaging.com
	Symbol for «Manufacturing year»		
	Complies with European directive 93/42/EEC		
	Symbol for «Type BF device»		
	Electronic and electrical equipment put on the market after 13/08/2005. This symbol indicates that this product should not be treated with household waste.		
	Symbol for «Equipotential plug»		
	Symbol for «UL/CSA timed fuses»		
	Symbol for «Read the instructions for use»		
	Symbol for «HDMI video output»		
	Symbol for «USB output»		

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