**delmont**  
imaging

**EN - Instructions for use**  
Instrumentation

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# TABLE OF CONTENTS

|                                                                      |           |
|----------------------------------------------------------------------|-----------|
| <b>1. Device description.....</b>                                    | <b>3</b>  |
| 1.1. Intended use                                                    |           |
| 1.2. User training and experience                                    |           |
| 1.3. Please note the following before use                            |           |
| 1.4. Function check                                                  |           |
| <b>2. Use.....</b>                                                   | <b>3</b>  |
| 2.1. Introducing the instrument into the endoscope                   |           |
| 2.2. Removing the tissues sample or foreign body                     |           |
| 2.3. Cutting tissue or sutures                                       |           |
| 2.4. Rescuing tissue or foreign bodies using baskets                 |           |
| 2.5. Withdrawing the instrument from the endoscope                   |           |
| <b>3. Preparation: Cleaning, disinfection and sterilization.....</b> | <b>5</b>  |
| 3.1. Warning and precautions                                         |           |
| 3.2. Cleaning and disinfection                                       |           |
| 3.3. Checking                                                        |           |
| 3.4. Maintenance                                                     |           |
| 3.5. Packaging                                                       |           |
| 3.6. Sterilization                                                   |           |
| 3.7. Storage                                                         |           |
| 3.8. Material resistance                                             |           |
| 3.9. Reusability                                                     |           |
| <b>4. Product type specific instructions.....</b>                    | <b>9</b>  |
| <b>5. Used symbols.....</b>                                          | <b>10</b> |

**REF** This manual is for references: D300 120 000 to D300 120 007.

## 1. Device description

### 1.1. Intended use

Biopsy forceps are used for the endoscopic removal of tissue for histological examination. Foreign body forceps and foreign body grabbers are used for the removal of foreign bodies. Scissors are used to cut tissue or sutures. Baskets are used to salvage tissue or foreign bodies.

These forceps have been developed exclusively for use in the gynaecology in combination with a suitable endoscope.

### 1.2. User training and experience

This instrument must only be used by a doctor or member of the clinical team under the supervision of a doctor with adequate training and experience in clinical endoscopic techniques. The user must have had sufficient training in endoscopic techniques. These instructions for use do not contain any description or explanation of endoscopic procedures.

### 1.3. Please note the following before use

Check the instrument for damage before every use. If the instrument is damaged or exhibits any other irregularities, the safety of the user or patient may be compromised. Injuries such as perforations, bleeding or injuries to the mucosa, as well as damage to the endoscope itself, may occur. Only use fully intact instruments. Always have a replacement instrument available.

 *This instrument has not been sterilised before shipping. It must therefore be prepared as described in Section 3 "Preparation: Cleaning, disinfection and sterilisation".*

### 1.4. Function check

 *If you discover a malfunction of or damage to the instrument of any kind, do not use the instrument. Use a replacement instrument instead. Damage or irregularities can compromise the safety of the patient or user, represent an infection risk, cause tissue irritation, perforation, bleeding, mucosal tissue injuries or serious damage to the equipment.*

 *Since infectious material from the patient represents an infection risk, appropriate protective clothing such as goggles, mouth masks, waterproof clothing and chemical-resistant gloves of an appropriate size and length to cover the skin must be used.*

## 2. Use

### 2.1. Introducing the instrument into the endoscope

Pull on the slider to close the jaws. Guide the instrument with the jaws closed into the endoscope's biopsy valve.

 *When introducing it into the endoscope, hold the instrument in as much the same axis as the endoscope as possible. Otherwise damage can occur to the instrument and/or endoscope.*

 Ensure that the instrument being used fits in the endoscope's working channel. Grip the instrument approximately 3 cm up from the working channel and advance it slowly and carefully with the jaws closed until the instrument clearly appears in the viewing field. Take care when the instrument exits the endoscope, since injuries such as perforations, bleeding and injuries to the mucosal tissue can be inflicted on the patient or the endoscope can be damaged.

 Hold the slider firmly when introducing the instrument into the endoscope. Otherwise the instrument can open and suddenly exit the endoscope. This can cause injuries to the patient such as perforations, bleeding or mucosal tissue injuries. The endoscope and / or instrument can also be damaged.

 Do not insert the instrument against resistance. Reduce the angular deflection until the instrument can be inserted effortlessly. Failing to do so can cause injuries to the patient such as perforations, bleeding or mucosal tissue injuries. The endoscope and / or instrument can also be damaged.

 Never advance the instrument suddenly. Doing so can cause injuries to the patient such as perforations, bleeding or mucosal tissue injuries. The endoscope and / or instrument can also be damaged.

## 2.2. Removing the tissue sample or foreign body

Guide the instrument to the site from which the tissue sample or foreign body is to be removed. Open the jaws and grasp the tissue or object to be removed.

 Do not push too hard, since doing so can cause injuries to the patient such as perforations, bleeding or mucosal tissue injuries.

## 2.3. Cutting tissue or sutures

Guide the instrument to the site at which the tissue or the sutures are to be cut. Open and close the instrument to separate the tissue or the sutures.

 Do not push too hard, since doing so can cause injuries to the patient such as perforations, bleeding or mucosal tissue injuries.

## 2.4. Rescuing tissue or foreign bodies using baskets

Guide the basket to the site from which the tissue or foreign body is to be removed. Open the basket, grasp the tissue or foreign body to be removed and carefully close the basket.

 Do not push too hard, since doing so can cause injuries to the patient such as perforations, bleeding or mucosal tissue injuries.

## 2.5. Withdrawing the instrument from the endoscope

 Actuate the handle to close the instrument. Slowly and carefully withdraw the instrument in a closed state out of the endoscope. Otherwise blood, mucus and other bodily fluids could escape and represent an infection risk.

 Should you feel resistance, reduce the angular deflection of the endoscope until the instrument can be withdrawn effortlessly. Forcible withdrawal can damage the instrument and / or endoscope. If the instrument no longer close fully, remove the endoscope with the instrument.

 *Ensure that no infectious material is lost as this can represent an infection risk for the patient and operator.*

## 3. Preparation: Cleaning, disinfection and sterilization

### 3.1. Warning and precautions

All instruments must be cleaned, disinfected and sterilised before every use. This is especially also true for the first use after delivery, since all instruments are shipped non-sterile: clean and disinfect after removal of the protective transport packaging, then sterilise after packaging. Effective cleaning and disinfection are an essential requirement for effective sterilisation.

The reprocessing of endoscopic equipment must remain the special responsibility of one person or several people. These individuals must be thoroughly trained in reprocessing methods and be aware of and fully understand the following points:

- your hospital's reprocessing guidelines,
- health and safety guidelines as defined by employment law,
- national and local hospital guidelines and principles ,
- the instructions in these instructions for use,
- the mechanical aspects of the instrument,
- the guidelines on hospital hygiene set down by the Robert Koch Institute (Germany),
- additional national specifications that must be complied with in relation to effective prion inactivation (e.g. mandatory use of a highly alkaline cleaning agent, prolonged sterilisation times).

Please note as part of your responsibility for the sterility of instruments in use:

- that essentially only methods for cleaning and disinfection and sterilisation are used that are adequately validated for the specific equipment and product,
- that the equipment used (CDD, steriliser) is regularly maintained and checked,
- that the validated parameters are maintained during each cycle.

Please note any other applicable legal requirements in your country as well as the hygiene regulations in force in your clinical practice or hospital. This is especially true for the various specifications governing effective prion inactivation (not applicable in the USA).

 *During use, the instrument comes into contact with potentially infectious tissue. To minimise the risk of infection, the instrument must be carefully and properly cleaned and disinfected after each use to remove micro-organisms and organic material. If an instrument is not adequately cleaned and disinfected, effective sterilisation is not possible. During reprocessing, always wear suitable protective clothing such as goggles, mouth protection, waterproof protective clothing and chemical-resistant gloves, since you may be exposed to infectious material and potentially hazardous chemicals.*

The reprocessing of instruments must - within the deadlines set out below - be commenced immediately after use, since dried-on organic residue can hamper effective cleaning, disinfection and sterilisation.

### 3.2. Cleaning and disinfection

For cleaning and disinfection purposes, a mechanical method (CDD: cleaning and disinfection device) should be used wherever possible. A manual method - even using an ultrasound bath - should only be used if a mechanical method is not available due to the significantly inferior effectiveness and reproducibility. The use of a manual cleaning and disinfection method must be safeguarded by additional product and method-specific validation which is the responsibility of the user.

#### **Pre-treatment must be carried out in both cases:**

Coarse soiling must be removed from the instrument immediately after use (within a maximum of 2 hours). Depending on the product type, the specific instructions set out in Section 4.0 must be followed.

#### **Procedure:**

1. With the instrument immersed in cleaning solution, remove any soiling visible on the jaws using a clean soft brush (or a clean, soft and lint-free cloth) which is used solely for this purpose. Never use metal brushes or wire wool.
2. Open and close the jaws, immersed in cleaning solution, at least 3 times.
3. If applicable: dismantle the instruments as much as possible.
4. Place the instrument with the jaws open, along with the handle, in an ultrasound bath filled with cleaning solution (the duration and concentration used in this step should be taken from the cleaning solution manufacturer's specifications).
5. Rinse the instruments for at least 1 minute under running water (temperature < 35°C/95°F).

 *If - for example for workplace' safety reasons - you are using a cleaning and disinfectant agent for this, please be aware that it must be free of aldehydes (otherwise blood-based soiling may become fixed), must demonstrate tested effectiveness (e.g. VAH/DGHM or FDA/EPA approval / clearance / registration or bear the CE mark), be suitable for the disinfection of medical instruments and be compatible with the instrument (see section entitled "Material resistance"). Please note that disinfectant products used during pre-treatment are used solely to protect personnel and do not replace the later disinfection stage that must be carried out after cleaning.*

#### **Mechanical cleaning / disinfection (CDD (Cleaning and Disinfection Device)):**

When choosing a CDD, bear in mind:

- that the CDD must essentially have a tested effectiveness (e.g. DGHM or FDA approval / clearance / registration or bear the CE mark in accordance with DIN EN ISO 15883),
- that, where possible, a thermal disinfection programme should be used ( $A_0$  value > 3000 or, on older equipment, at least 5 min. at 90°C/194°F). With chemical disinfection there is a risk of disinfectant residue on the instruments,
- that the programme used is suitable for the instruments and includes adequate rinsing cycles,
- that only sterile or low-germ (max. 10 bacteria/ml) and low-endotoxin (max. 0.25 endotoxin units / ml) water is used (e.g. purified water / highly purified water),
- that the air used for drying is filtered (oil free, low-germ and low-particle),
- that the CDD is maintained and checked regularly.

When choosing the cleaning agent system to be used, bear in mind:

- that it must essentially be suitable for the cleaning of instruments made from metal and plastic,
- that - where no thermal disinfection method is used - a suitable disinfectant agent with tested effectiveness (e.g. VAH/DGHM or FDA/EPA approval / clearance / registration or bear the CE mark) is used and that this is compatible with cleaning materials being used,

- that the chemicals used are compatible with the instruments (see section entitled "Material resistance").

The concentrations, temperatures and exposure times specified by the cleaning and possibly disinfectant agent manufacturer, as well as rinsing specifications, must be complied with at all times.

Depending on the product type, the specific instructions set out in Section 4.0 must be followed.

## Procedure:

1. Place the dismantled instruments with the jaws open in the CDD. Ensure that the instruments are not touching.
2. If applicable: Place the handle parts in a small parts basket in the CDD.
3. Start the programme.
4. Remove the instruments from the CDD when the programme finishes.
5. Check and package the instruments as soon as possible after removal (see sections entitled "Checking", "Maintenance" and "Packaging", if necessary after additional drying in a suitable location).

 *Evidence of the fundamental suitability of the instruments for effective mechanical cleaning and disinfection has been provided by an independent, accredited testing laboratory using the CDD G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the cleaning agent Neodisher medizym (Dr. Weigert GmbH & Co. KG, Hamburg). The method described above was applied in this context.*

## 3.3. Checking

Check the instrument's function by opening and closing the jaws several times. Check all of the instruments after cleaning or cleaning / disinfection for corrosion, damaged surfaces, cracks, soiling and discolouration and remove any damaged instruments from use (the number of times an instrument can be reused is limited - see the section entitled "Reusability"). Instruments that are still soiled must be cleaned and disinfected again.

## 3.4. Maintenance

Instrument oils and greases must not be used.

## 3.5. Packaging

Place the cleaned and disinfected instruments coiled up in a single-use sterilisation pack (single or double packaging) and / or sterilisation containers. The specific instructions set out in Section 4.0 must be followed.

The packaging must meet the following requirements (material / process):

- DIN EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance),
- suitable for steam sterilisation (temperature resistance up to at least 138°C (280 °F) with adequate steam permeability),
- adequate protection of the instruments or sterilisation packaging against mechanical damage,
- regular maintenance in accordance with the manufacturer's specifications (sterilisation containers).

Please also note the instructions for use on sterile packaging and the sealing device.

 *Ensure that the sterile packaging is not damaged when the instrument is inserted into it. Do not bend rigid tube shaft instruments. Do not use excessive force and ensure that the shaft is not bent.*

### 3.6. Sterilization

Only the sterilization methods listed below must be used for sterilization. Other sterilization methods are not permitted.

The specific instructions set out in Section 4.0 must be followed.

#### **Steam sterilization:**

- Fractionated vacuum method: at least three vacuum steps with adequate product drying: the actual drying time required depends directly on the parameters that are the sole responsibility of the user (loading configuration and density, steriliser condition, etc.) and must therefore be determined by the user. Regardless of this, drying times must not be less than 20 minutes.

*The use of the less-effective gravitation method is not permitted.*

- Steam steriliser compliant with DIN EN 13060 / DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance).
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ (picking) and product-specific performance assessment (PQ)).
- Maximum sterilization temperature 134°C (273°F; plus tolerance as defined in DIN EN ISO 17665).
- Sterilization time (exposure time at sterilisation temperature) of at least 20min. at 121°C (250°F) or at least 5min. (or 18min for prion inactivation, not relevant for the USA) at 132°C (270°F)/134°C (273°F).

*Evidence of the fundamental suitability of the products for effective steam sterilization has been provided by an independent, accredited testing laboratory using steam steriliser HST 6x6x6 (Zirbus technology GmbH, Bad Grund) and the fractionated vacuum method. Typical conditions found in hospitals and medical practices, as well as the method described above, were used.*

The flash sterilization method is not permitted.

Do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization or plasma sterilization.

*After sterilization, check the sterile packaging for damage, adequate drying and to ensure that the seal is intact.*

*If a defect is found, the packaging and sterilization steps must be repeated.*

### 3.7. Storage

After sterilization, the instruments must be stored in the sterilization packaging in a dry and dust-free environment. Keep the instrument at room temperature in a clean, dry place and protect it from direct sunlight.

Please also note the instructions given by the sterile packaging manufacturer.

### 3.8. Material resistance

When choosing cleaning and disinfection agents, please ensure that they do not contain the following:

*(Please add material data sheets and user reports as appropriate, possible chemical substances to be avoided in standard products)*

- strong organic, mineral or oxidising acids (lowest permissible pH value 5.5),
- strong alkalis (highest permissible pH 11, neutral / enzymatic or slightly alkaline cleaners are recommended),
- organic solvents (e.g. ether, ketones, benzines), fluorinated alcohols,

- oxidising agents (e.g. hydrogen peroxide),
- halogens (chlorine, iodine, bromine),
- aromatic / halogenated hydrocarbons,
- formamide,
- trichloroethylene / perchloroethylene.

Never clean any instruments with metal brushes or wire wool.

Do not expose any instrument to temperatures higher than 138°C (280°F).

### 3.9. Reusability

The instruments can - with appropriate care and provided they are undamaged and unsoiled - be reused up to 100 times; any further use beyond this or the use of damaged and / or soiled instruments is the responsibility of the user.

## 4. Product type specific instructions

| Handle                                                                              | In addition to Section 3.2. Pre-treatment                                                                                                                                                                                                                                                                                                                      | In addition to Section 3.2. Mechanical cleaning / disinfection                                | In addition to Section 3.5 Packaging                                                                        | In addition to Section 3.6 Sterilisation                                                                    |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|
|  | <p>3. Dismantle the forceps top part from the handle as follows:</p> <ul style="list-style-type: none"> <li>- Unscrew the locking screw</li> <li>- Unhook the handle section from the forceps top part</li> <li>- Unscrew the screw sleeve on the handle</li> <li>- Pull the forceps top part out of the handle</li> </ul> <p>Otherwise standard procedure</p> | <p>2. Handle parts in a small parts basket in the CDD</p> <p>Otherwise standard procedure</p> | <p>Attach the handle, but with the screw sleeve on the handle loose</p> <p>Otherwise standard procedure</p> | <p>Attach the handle, but with the screw sleeve on the handle loose</p> <p>Otherwise standard procedure</p> |

No liability is accepted for failure to comply with this procedure.

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