



**delmont**  
imaging


**Our colposcopes**  
now combined with **iCare<sup>mini</sup>**  
for even greater performance

# 1. THE PRINCIPLE

---


# What is colposcopy?

- **Colposcopy** is an in-depth examination that allows **direct visualization** of:
  - The cervix
  - The vagina walls
- It is carried out using a **binocular magnifying instrument (microscope)** combined with a light source.
- The aim of this examination is to localize **any abnormal lesions**, and consequently to locate **sites for biopsies**.
- Colposcopy is a relatively painless examination that requires **no anaesthetic**.



# HPV, the human papillomavirus

- The vast majority of **cervical cancers** are caused by the **human papillomavirus (HPV)**, which is sexually transmitted. It is a very common infection among sexually active men and women.
- There are over **100 different types** of HPV:
  - In most cases, HPV causes no health problems and disappears naturally.
  - Some types only cause genital warts.
  - But sometimes a woman is infected with a "**high-risk**" HPV (**~15 types**) that doesn't go away on its own.
  - And in rare cases, an HPV infection that doesn't disappear **turns into cervical cancer**.
- Cervical cancer takes **10 to 20 years to develop** following HPV infection.



# Colposcopy in screening


- The **process of screening for uterine cancer** and the recommendations of medical societies vary from country to country.
- Overall, however, **colposcopy is recommended** in the following cases:

- A positive HPV test
- One pathological PAP (papanicolaou) test
- Two consecutive unsatisfactory PAP tests.
- Follow-up examination after surgery.

An **HPV test** uses molecular research to detect the presence or absence of **viral DNA**.

A **PAP test** is the microscopic examination (cytology) of the sample **for abnormal or precancerous cells**.


## Cervico-uterine smear test




- Cervical cancer is one of the most common cancers in women. But it's also **one of the easiest to prevent!**

# The different steps of the examination

- After inserting a speculum to keep the vaginal walls apart, and positioning the colposcope, the doctor follows the steps below:
  1. Visualization of the cervix without dye or filter for an initial **general assessment** (secretions, inflammation, outgrowth, etc.).
  2. Use of the **green filter** for **vascularization analysis**. The cervix must be cleaned with saline beforehand.
  3. **Acetic acid test**: vinegar in very dilute concentration which whitens tissues with a high protein load, i.e. abnormal cells called "acidophilic".
  4. **Lugol test**: iodine base that blackens normal mucous membranes. Complementary to acetic acid, the test confirms a diagnosis and better defines the limits of abnormal areas.
  5. Depending on the results, a **cervical biopsy and/or endocervical curettage** may be performed.
  6. Examination of the **vaginal walls** (the acetic acid test can also be performed on these areas).



No dyes or filters



Green filter



Acetic acid test



Lugol test

# Different locations and types of use

- Depending on the purpose of the colposcopy (diagnostic or operative), the examination will not be performed under the same conditions:



## DIAGNOSTIC COLPOSCOPY

- In the consultation room, whether in private practice or in a healthcare center (private or public).
- No anesthesia.
- No advanced surgical interventions.
- Examination carried out with the primary aim of locating biopsy sites and performing them, as described above.



## OPERATIVE COLPOSCOPY

- Procedure performed **when HPV lesions are detected** following the results of a biopsy performed during a previous diagnostic colposcopy.
- In the operating room of a healthcare center (private or public).
- With general or locoregional anesthesia.
- Treatment options to remove lesions include cryotherapy, laser ablation, thermoelectric loop or scalpel conization.

# 2. RANGE OVERVIEW

---



# Three ergonomics available

## Vertical post colposcope



- Suitable for **small rooms**.
- The **least expensive** in the range.
- 3 magnification levels.

## Swing-arm colposcope



- High **comfort of use**.
- The most **popular**.
- 5 magnification levels.

## Over-the-shoulder colposcope




- Ergonomics designed for the **operating room**.
- 5 magnification levels

# Two types of lens available


## Fixed focus lens

- The colposcope must be positioned at the exact distance of the lens focal length (usually 30cm for colposcopy).
- **Focus range of 11mm.**
- The colposcope must be moved regularly during the examination.



## Variable-focus lens

- The colposcope can be positioned without high precision (ideally ~30cm from the cervix).
- **Focus range from 20 to 35cm.**
- The colposcope does not need to be moved during the examination.




# Optical customization possible

Three eyepiece models available:

- 10x / 18mm
- **12,5x / 16mm (by default)**
- 16x / 16mm

Three binocular models available:

- **45° binocular (by default)**
- 0° binocular
- Tiltable binocular from 0° to 220°





**Regarding the eyepieces:**

- The different models have an impact on **the diameter of the field of vision.**
- **Optional reticles can be added** inside them at the time of purchase.

Two types of lens available:

- **Fixed-focus** lens: 175mm, 200mm, 250mm, 300mm or 400mm.
- **Variable-focus** lens: Range from 200 to 350mm

# The beam splitter



- As soon as an additional imaging device is to be connected, a **beam splitter** is required.
- Our beam splitter is **double**: it divides each of the two initial optical beams in two, to enable the connection of two separate devices.
- Division is made possible by the use of **prisms**.
- Depending on the **prism's coating**, the incident light beam will not be equally separated between the two beams created.
- In our case, our beam splitter consists of:
  - a **50/50** splitting prism
  - an **80/20** splitting prism

# Connecting a DSLR camera



- To connect a DSLR camera to the colposcope, a **special camera adapter** must be attached to the beam splitter.
- **Three connector versions** are available:
  - Canon brand
  - Nikon brand
  - Sony brand



DSLR cameras require a significant amount of light to produce resolute photos. We recommend installing them on the **50/50 connector** of the beam splitter.

# Connecting a video camera



- To connect a video camera to the colposcope, a **specific video adapter** must be attached to the beam splitter (80/20 connector).
- **The connection of the adaptor is CS-mount**, but a ring is available to obtain a C-mount.

## Dedicated colposcope camera

- Direct connection to video adapter.
- Several types of video camera available:



## Camera not dedicated to colposcope

- For example, when a hysteroscopy column is already present in the room.
- **Add a bonnette adapter** for easier use.



# Adding a monitor



- A monitor support is available in the catalog for **integrated video display** on the colposcope:
  - Compatible monitor size: up to 24".
  - VESA mount 100x100
- An image capture device with video output must be connected to the beam splitter (DSLR camera, iCare mini, etc.).




This monitor support **is not compatible with the vertical post colposcope.**

# Micromanipulator

- A micromanipulator **couple a surgical CO<sub>2</sub> laser to the colposcope's vision** for the treatment of uterine lesions.
- The joystick moves the laser beam to the desired target.
- A **laser adapter** is also required, but will differ according to the brand and model of the customer's laser:

Reference	Designation
D100 300 511	Micromanipulator for Delmont colposcope
D100 300 013	CO <sub>2</sub> laser adapter for micromanipulator





# Armchair mounting



- The colposcope can be mounted on a chair frame using an **additional mounting arm**.
- Only **Schmitz gynecological chairs** are compatible with our mounting arm.
- This configuration **considerably reduces overall footprint**.



**Only swing-arm colposcope ergonomics are compatible** with this installation.

# And even more accessories

## Protective cover

- For use in **operating rooms**.
- **Sterile**, single-use device.
- Allows the user to manipulate the colposcope head while wearing sterile clothing.



## Co-observation arm

- Allows **two separate operators to look through** the colposcope simultaneously.
- Two versions are available:



**2D assistant arm**





**3D assistant arm**

# 3. STRENGTHS

---

# High-definition optical quality

- The primary aim of colposcopy is to **make a correct diagnosis, and to detect potential uterine cancer as early as possible.**
- Good quality of vision is therefore the **most important feature.**
- Our colposcopes incorporate **very high-definition optical components.**
- And our lenses are **apochromatic** (chromatic aberration correction).



# Complete and ultra-customizable range

- Three different ergonomics, a wide range of accessories: **we have the perfect solution for every customer.**
- **Upgradeable at any time and for any accessory!**



We launched the range in 2020. **In France alone**, by the end of 2022 we had sold **57 colposcopes**. And in 2023 alone, we have sold **42 units**.  
The ratio between the three ergonomics is as follows:



57%



38%



5%

# Variofocus, optimised comfort



- The variable-focus lens **significantly increases the gynecologist's ease of use.**
- Having a **sharp image throughout the examination** becomes extremely easy.
- Ultimately, this ease of use translates into **significant time savings** for the doctor, for each examination.




We sell **89%** of our colposcopes with a variofocus lens.

# Colposcopy report

- Via our imagyn application, available on iPad, a **colposcopy report** is available with many advantages:
  - Integration of exam photos.
  - Anatomical drawings that you can annotate.
  - Complete list of fields, fillable in just a few clicks.



To use the imagyn application, customers must also **have iCare or iCare mini.**



# 4. THE USE

---



# Head operation

**Magnification knob**  
available on either side  
(3 or 5 magnifications  
depending on model)

**Beam splitter** for  
connecting accessories

**Integrated mechanical  
green filter.**  
Push/Pull tab.

**Handle for positioning  
the colposcope head**



**Interpupillary distance and  
diopter adjustment**

**Focus knob** 360° positionable:


- Variofocus: range from 200 to 350mm
- Fixed lens: 11mm range

Setting the **mobility hardness**  
of the colposcope head


**Light cable** connection

# Colposcope body operation

## Vertical post colposcope



## Swing-arm / Over-the-shoulder colposcope



# How to set up a colposcope?

1.

**Turn ON the light source**  
Switch ON the main button.  
Once the light source is ready, press the Lamp button.



2.

**Inter-pupil space setup**  
Look through both eyepieces.  
Adjust the distance manually until each eye image overlaps.



3.

**Diopter adjustment (1/4)**  
**Only for myopia and/or hyperopia.**  
Set each eyepiece to 0 (line to line). Target the view on an object.



4.

**Diopter adjustment (2/4)**  
Set the magnification to 5.  
Make a rough focus by moving the head and a precise one with the knob.



5.

**Diopter adjustment (3/4)**  
Without moving the colposcope position or the focus knob, set the magnification to 1.




6.

**Diopter adjustment (4/4)**  
Adjust both eyepiece diopters by turning them until you reach a sharp image.  
**Record your values for future use.**



# Use with a video camera

- Set up the colposcope as shown in the previous slide.
- Target a white surface with the colposcope head and **launch a white balance** via the camera head.
- Place the colposcope ~30 cm from the cervix and **focus with the objective knob (1)** while looking through the eyepieces.
- Then look at the image on the monitor and **focus using the video adapter ring (2)**.
- If the bonnette adapter (3) is installed, **focus using the two rings (2) and (4)** to obtain sharpness on the screen.
- For the rest of the examination, **use only the lens knob (1)**. Rings (2) and (4) must not be used!



# iCare mini buttons

- The iCare mini buttons allow you to:

- Short press: **take photos**
- Long press: start/stop **video recording**

Launch **white balance** (long press)



- Short press: activates or deactivates **image freeze**
- Long press: launch or exit **preview mode**

Video recording status indicator (**flashes red** when recording)

Switch between **standby and exam mode**


## Product status indicator light :

- Rapid white flash: booting phase (~30sec)
- Slow white flash: standby mode
- Steady white light: examination mode



**The three main buttons are also used to navigate** the home and settings screens.

# iCare mini connectors



**Wi-Fi antennae built into the device** for communication with imagyn: transfer of streaming video and photos of the examination.

Three connectors available on the underside:

- 1 USB 2.0 output for connecting a **USB storage key for photos and videos.**
- 1 USB 2.0 output for connecting a **Wi-Fi dongle** to connect the camera to the healthcare centre network.
- **1 micro-HDMI video output** for connection to a monitor.

On the left-hand side, there is a **Jack input for connection to an external power supply.**  
*Our colposcopes incorporate an external power supply in the LED source housing.*




# iCare mini interface

## Home and settings screen




- Several settings are possible:
  - **iPad connection configuration**
  - **Language:** 14 languages available
  - **Image:** choice of 3 presets or manual adjustment
  - **USB key:** current capacity and deleting option
  - **General information:** launch camera software updates and factory reset.

## Features in Exam mode



➤ Taking photos

➤ Video recording



➤ Photo preview mode

# 5. THE COMPETITION

---





KP 3000




SOM 52



- **Motorized** magnification and focus available on the swing-arm version (SOM 52).

- No **variable-focus lens** available (variofocus).
- **Overall inferior quality:** weight and stability of colposcope stand, quality of optical components, etc.
- **Unintuitive software**, with limited reporting quality.






iView 21 or 31



- **Specific H.A.S.I filter** for better contrast than conventional green filters.
- **Measurement scale overlay** available at any time.

- **Limited range of ergonomics available:** stand-mounted or chair-mounted swing-arm versions.
- **No upgradability of image capture possible with Model 21.**
- **Unintuitive software,** with limited reporting quality.





CooperSurgical®

**LEISEGANG**  
The Colposcope Company

OptiK Model 1



- **Reticles available** by default on all their models for lesion measurement.

- No **variable-focus lens** available (variofocus).
- "**Old-generation**" design.
- No **model upgradability**.
- All models **only have a maximum of three magnifications levels**.
- **Overall inferior quality:** weight and stability of colposcope stand, quality of optical components, etc.
- **Panel PC software** with limited reporting quality (captioned photos and videos only).



OptiK Model 2



OP-C2L



OP-C2




OP-C5




- **Upgradability of image capture possible on all models.**

- **Only one ergonomic design available: swing-arm colposcope.**
- **Variable-focus lens available only on the most expensive model.**
- **Overall inferior quality: weight and stability of colposcope stand, quality of optical components, etc.**
- **Old-generation software, not very intuitive, with limited reporting quality.**




# Centrel

The Clear Vision



- **Continuous magnification available only on the most expensive model (Z4).**

- **No variable-focus lens available (variofocus).**
- **Overall inferior quality:** weight and stability of colposcope stand, quality of optical components, etc.
- **Old-generation software,** not very intuitive, with limited reporting quality.




Kv 20

ISIDE BETA

ISIDE.hub

- A **wide choice** of models with different image capture solutions available.

- **No variable-focus lens** available (variofocus).
- **No model upgradability.**
- Still not very popular in gynecology, and highly variable from country to country (Italian manufacturer). A company **initially focused on dermatology.**
- Software on PC or Panel PC (depending on model) with **limited reporting quality.**



# Video-colposcopes

- The video-colposcopy market is a completely different market (2D), for very specific users. **It should not be considered as direct competition.**
- Some examples of brands and models:

Schmitz – Vidan 2



KAPS – ViCo S HD



Duomed – Dysis



Edan – C3A




# 6. AND FINALLY...


---




# Documents to have




Sales leaflet




Example of a report



Datasheet



Patient information sheet



Quick-start sheet



All documents are available on our website and Extranet.