

LIST OF DISINFECTANTS FOR ENDOSCOPES



Cleaning/Disinfection	Detergent	Laboratory
Disinfection	Actanios HLD	Anios
Disinfection	Cides oPA	Johnson& Johnson
Disinfection	Disinfectand	Olympus
Disinfection	Gigasept AF forte	Schülke
Disinfection	Gigasept FF (neu)	Schülke
Disinfection	Gigasept Instru AF	Schülke
Disinfection	Kombi Instru N	Antiseptica
Disinfection	Korsolex Extra	Bode
Disinfection	Korsolex Extra	Bode
Disinfection	Micro-Cide 28	Micro-Scientific
Disinfection	Mikrozid AF liquid	Schülke
Disinfection	Neodisher endo DIS active	Dr. Weigert
Disinfection	Neodisher Endosept GA	Dr. Weigert
Disinfection	Neodisher Mediklar	Dr. Weigert
Disinfection	Neodisher Septo DA	Dr. Weigert
Disinfection	Neodisher Septo DN	Dr. Weigert
Disinfection	Neodisher Septoclean	Dr. Weigert
Disinfection	Rely+On Perasafe	Lanxess
Disinfection	Sekusept aktiv	Ecolab
Disinfection	Thermosept DK	Schülke
Disinfection	Thermosept ED	Schülke
Disinfection	Thermosept ER	Schülke
Disinfection	THP protect	Schülke
Disinfection	Triacid N	Antiseptica
Disinfection	Tristel Wipes	Tristel
Cleaning	ProEZ foam	Certol
Cleaning	Actanios LDI	Anios
Cleaning	Bodedex forte	Bode
Cleaning	Cleaner	Olympus
Cleaning	Deko Neutral Detergent	Dekomed
Cleaning	Deko Rinse	Dekomed
Cleaning	Endozime Xtreme Power	Ruhof

LIST OF DISINFECTANTS FOR ENDOSCOPES



COMPATIBLE

Cleaning/Disinfection	Detergent	Laboratory
Cleaning	Neodisher Endoclean	Dr. Weigert
Cleaning	Neodisher Mediclean	Dr. Weigert
Cleaning	Neodisher Mediclean forte	Dr. Weigert
Cleaning	Neodisher Mediklar	Dr. Weigert
Cleaning	Prolystica® 2X Concentrate	Steris
Cleaning	Sekusept Multienzyme	Ecolab
Cleaning	Sekuseptforte S	Ecolab
Cleaning	Sprayclean 2000	Olympic Chemical
Cleaning/ Disinfection	neodisher Septo Active	Dr. Weigert
Sterilization	Autoklav 134 °C/ 7 Min	N/A
Sterilization	Sterrad	Johnson& Johnson



NON COMPATIBLE

Cleaning/Disinfection	Detergent	Laboratory
Disinfection	Trioton Rapid AF	Pliwa
Disinfection	Alka 100	Alkapharm
Disinfection	Cleanisept Wipes	Dr. Schuhmacher
Disinfection	excel D	Anios
Disinfection	Hexanios G+R	Anios
Disinfection	HMI Peroster CE	HMI
Disinfection	REDA Multi-zyme	REDA
Disinfection	Sekusept plus	Ecolab
Disinfection	Steranios 2%	Anios
Cleaning	Instrunet 4EZ+T	Inibsa
Sterilization	Reno-S30	Renosem