

HF cables



Photo	Reference	Description
Bipolar resection cable (ReseCare & MaxiCare)		
	D200 190 000	Covidien bipolar resection cable. Compatible with the Valleylab ForceTriad generator.
	D200 190 001	Erbe bipolar resection cable. Compatible with Erbe Vio 300 and Storz Autocon III 400 generators.
	D200 190 008	Bowa bipolar resection cable. Compatible with Bowa Arc 350 and 400 and Storz Autocon III 400 generators.
	D200 190 006	2 pins - 29mm bipolar resection cable. Compatible with all 2 pins - 29mm connection generators such as Erbe ES350, Pulsar MB 350, Megadyne and Lamidey Noury.
	D200 190 007	Emed (SDS plug) bipolar resection cable. Compatible with the Emed Spectrum generator.
	D200 190 012	Berchtold or Martin bipolar resection cable. Compatible with Berchtold, Aesculap and Martin Maxium generators.

Photo	Reference	Description
Monopolar resection cable (ReseCare)		
	D200 190 005	Valleylab monopolar resection cable. Compatible with Bovie, Valleylab and Conmed generators.
	D200 190 004	Erbe T-series monopolar resection cable. Compatible with Erbe T-serie, Martin and Berchtold generators (banana plug).
	D200 190 011	Berchtold monopolar resection cable. Compatible with Martin, Berchtold, Aesculap and Rudolf generators.
	D200 190 010	Erbe ACC/ICC monopolar resection cable. Compatible with Erbe ACC/ICC and Storz generators.
HF cable for bipolar electrodes (EasyCare)		
	D200 190 003	Valleylab cable for bipolar electrode. Compatible with Valleylab generators and international models.
	D200 190 002	Erbe cable for bipolar electrode. Compatible with Erbe and Storz generators.
	D200 190 009	Berchtold cable for bipolar electrode. Compatible with Berchtold and Martin generators.