

Product portfolio 2019

1. Equipment

Complete solutions		
Photo	Reference	Description
	D400 100 501	iCare+ . Imaging station for hysteroscopy, including: - iCare (coupler included), - iPad Pro 12,9", 256Gb (including keyboard and stylus), - imagyn license, - Neovo 22" monitor, - Cart, - Cloud storage of 5 Gb for 1 year.
	D400 100 503	iCare+ . Same as D400 100 501 but including isolation transformer on the cart.

Cameras, couplers & zooms		
Photo	Reference	Description
	D400 100 500	iCare. HD camera including F22 coupler. Light cable connector of Storz type.
	D400 100 504	iCare. HD camera including F22 coupler. Light cable connector of Olympus type.
	D100 100 000	iCare camera control unit only. Light cable connector of Storz type.
	D100 100 001	iCare camera control unit only. Light cable connector of Olympus type.

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat


Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

Photo	Reference	Description
	D400 100 506	iCare camera head only, with F22 coupler.
	D100 110 000	iCare camera head only, without coupler.
	D100 120 000	Coupler of focal length 22mm.
	D100 120 001	Coupler of focal length 18mm.
	D100 120 002	Zoom coupler of focal length 16 to 34mm
	D900 000 119	iCare 16Gb USB stick.

Digital support, application & Cloud		
Photo	Reference	Description
	D100 200 501	iPad Pro 12.9", white, 256Gb, supplied with Apple pencil, keyboard and imagyn license
	D100 200 000	iPad Pro 12.9", white, 256Gb.
	D100 200 001	Apple pencil for iPad Pro.
	D100 200 002	Smart Keyboard for iPad Pro 12.9" - French (AZERTY)

Photo	Reference	Description
 imagyn	D100 200 101	imagyn license.
	D100 200 103	imagyn sub-license.
 iZap Control	D100 200 102	iZap Control license (to directly download on the Apple Store)
	D100 200 201	5 Gb in the medical Cloud for 1 year.
	D100 200 200	10 Gb in the medical Cloud for 1 year.
	D100 200 502	iCare - iPad wired connection kit.

Monitors		
Photo	Reference	Description
	D200 110 000	Neovo 22" monitor with optical glass. Full HD resolution.
	D200 110 001	Neovo 22" monitor, medical grade, with optical glass. Full HD resolution.

Carts		
Photo	Reference	Description
	D200 100 000	Cart. 2 fixed shelves, 1 infusion stand, 1 camera holder.



delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat


Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

Photo	Reference	Description
	D200 100 003	Cart with isolation transformer. 2 fixed shelves, 1 infusion stand, 1 camera holder.
	D200 100 500	Arm and stand for iPad Pro.
	D200 100 501	Arm and stand for iPad Pro. Secure version.
	D200 100 001	Double pivot arm.
	D100 200 003	Stand for iPad Pro.
	D100 200 004	Stand for iPad Pro. Secure version.

Irrigation pumps		
Photo	Reference	Description
	D110 100 011	iFlow. Irrigation pump. Without any function pre-installed
	D110 100 500	iFlow. Irrigation pump. Hysteroscopy function.
	D110 100 501	iFlow. Irrigation pump. Hysteroscopy and Laparoscopy function.

Photo	Reference	Description
	D110 100 004	Transponder for Laparoscopy function.
	D110 100 007	Transponder for Hysteroscopy function.
	D110 100 008	Transponder for Arthroscopy function.
	D110 100 009	Transponder for Urology function.
	D110 100 010	Double pedal foot switch.
	D110 100 000	Digital bottle scale for input/output balance.
	D110 100 006	Bemis containers for secretions. 3L. Delivered by 2.
	D110 100 001	Box of 10 sterile irrigation tubes. Single use.
	D110 100 002	Box of 10 sterile suction tubes - 2 connectors. Single use.
	D110 100 005	Irrigation tube reusable 20 times.
	D110 100 003	Suction vacuum tube with filter, non sterile, for 30 days of use.

Ultracision			
Photo	Reference	Description	
	D300 300 000	Ultrasound generator.	
	D300 300 001	Hand piece for ultrasonic shears.	
	D300 300 002	L: 36cm	iScalpel. 5,5mm ultrasonic shears. Sterile and single use. Delivered by 3 only.
	D300 300 003	L: 23cm	
	D300 300 004	L: 14cm	
	D300 300 005	Foot switch (MIN / MAX) to connect to the ultrasound generator.	

2. Instrumentation





Complete solutions		
Photo	Reference	Description
	D400 100 017	EasyCare. Complete hysteroscopy set, including: - HD hysteroscope 2.9mm 30°, L: 300mm autoclavable, - diagnostic sheath 4mm rotatable stopcock, - inner operating sheath 4.9mm, 7Fr. operative canal, - outer operating sheath 5.5mm rotatable stopcock, - graduated scissors 7Fr. semi rigid, L:340mm, single action jaws, with flushing port, - graduated forceps (alligator) 7Fr., semi rigid, L:340mm, double action jaws, with flushing port, - graduated biopsy forceps 7Fr., semi rigid, L:340mm, double action jaws, with flushing port, - graduated IUD grasping forceps 7Fr., semi rigid, L:340mm, double action jaws, with flushing port, - light cable, grey, Ø: 3.5mm, L: 2.3m, delivered with Storz connectors.
	D400 100 000	EasyCare. Instruments without flushing port version.
	D400 100 011	ReseCare. 18.5Fr. bipolar resection set, including: - HD hysteroscope 2.9mm, 12°, L:300mm, autoclavable, - passive working element for 18.5Fr. resectoscope, - outer sheath, 18.5Fr., double flow, - inner sheath, rotatable, 17.5Fr., - cutting loop 30°, - electrode rollerball, - electrode knife, - cutting loop, straight, - light cable, grey, Ø: 3.5mm, L: 2.3m, delivered with Storz connectors.
	D400 100 013	ReseCare. Active version & 12° endoscope.

Photo	Reference	Description
	 D400 100 001	ReseCare. Passive version & 30° endoscope.
	 D400 100 003	ReseCare. Active version & 30° endoscope.
	 D400 100 008	MaxiCare. 26Fr. bipolar resection set, including: - HD hysteroscope 4mm, 12°, L:300mm, autoclavable, - passive working element for 26Fr. resectoscope, - outer sheath, 26Fr., double flow, - inner sheath, rotatable, 24Fr., - cutting loop 30° - cutting loop straight - electrode rollerball - electrode knife - light cable, grey, Ø: 3.5mm, L: 2.3m, delivered with Storz connectors.
	 D400 100 004	MaxiCare. Active version & 12° endoscope.
	 D400 100 002	MaxiCare. Passive version & 30° endoscope.
	 D400 100 014	MaxiCare. Active version & 30° endoscope.



Endoscopes		
Photo	Reference	Description
Hysteroscopes		
	D300 100 003	HD hysteroscope 2.9mm 0° 300mm autoclavable.
	D300 100 004	HD hysteroscope 2.9mm 12° 300mm autoclavable.
	D300 100 000	HD hysteroscope 2.9mm 30° 300mm autoclavable.
	D300 100 005	HD hysteroscope 4mm 0° 300mm autoclavable.
	D300 100 002	HD hysteroscope 4mm 12° 300mm autoclavable.
	D300 100 001	HD hysteroscope 4mm 30° 300mm autoclavable.
	D300 100 007	HD hysteroscope 4mm 70° 300mm autoclavable.
Laparoscopes		
	D300 100 015	HD Mini laparoscope 2.9mm 0° 207mm autoclavable.
	D300 100 016	HD Mini laparoscope 2.9mm 30° 207mm autoclavable.
	D300 100 020	HD laparoscope 5.5mm 0° 300mm autoclavable HM.
	D300 100 021	HD laparoscope 5.5mm 30° 300mm autoclavable HM.
	D300 100 022	HD laparoscope 5.5mm 45° 300mm autoclavable HM.
	D300 100 026	HD bariatric laparoscope 5.5mm 0° 424mm autoclavable HM.
	D300 100 027	HD bariatric laparoscope 5.5mm 30° 424mm autoclavable HM.
	D300 100 028	HD bariatric laparoscope 5.5mm 45° 424mm autoclavable HM.

Photo	Reference	Description
	D300 100 032	HD laparoscope 10mm 0° 330mm autoclavable.
	D300 100 033	HD laparoscope 10mm 30° 330mm autoclavable.
	D300 100 034	HD laparoscope 10mm 45° 330mm autoclavable.
	D300 100 038	HD bariatric laparoscope 10mm 0° 450mm autoclavable.
	D300 100 039	HD bariatric laparoscope 10mm 30° 450mm autoclavable.
	D300 100 040	HD bariatric laparoscope 10mm 45° 450mm autoclavable.
	D300 100 041	Operating laparoscope 10mm, 0°, 270mm, non-autoclavable, 45° angulated eyepiece, with 5.1mm working channel.
	D300 100 043	Operating laparoscope 11mm, 0°, 270mm, non-autoclavable, 45° angulated eyepiece, with 6.1mm working channel.
	D300 100 042	Operating laparoscope 10mm, 0°, 270mm, non-autoclavable, 90° angulated eyepiece, with 5.1mm working channel.
	D300 100 044	Operating laparoscope 11mm, 0°, 270mm, non-autoclavable, 90° angulated eyepiece, with 6.1mm working channel.
Sinusscopes		
	D300 100 052	Sinuscope 2.7mm 0° 140mm autoclavable.
	D300 100 053	Sinuscope 2.7mm 30° 140mm autoclavable.
	D300 100 054	Sinuscope 2.7mm 70° 140mm autoclavable.
	D300 100 055	Sinuscope 2.7mm 0° 195mm autoclavable.
	D300 100 056	Sinuscope 2.7mm 30° 195mm autoclavable.
	D300 100 057	Sinuscope 2.7mm 70° 195mm autoclavable.

Photo	Reference	Description
	D300 100 058	Sinuscope 3mm 0° 150mm autoclavable.
	D300 100 059	Sinuscope 3mm 30° 150mm autoclavable.
	D300 100 060	Sinuscope 3mm 45° 150mm autoclavable.
	D300 100 061	Sinuscope 3mm 70° 150mm autoclavable.
	D300 100 062	HD sinuscope 4mm 0° 185mm autoclavable.
	D300 100 063	HD sinuscope 4mm 30° 185mm autoclavable.
	D300 100 064	HD sinuscope 4mm 45° 185mm autoclavable.
	D300 100 065	HD sinuscope 4mm 70° 185mm autoclavable.
Laryngoscopes		
	D300 100 073	Micro laryngoscope 4mm 70° 185mm autoclavable, reversed light post.
	D300 100 074	Laryngoscope 6mm 70° 160mm autoclavable HM.
	D300 100 075	Laryngoscope 6mm 70° 160mm autoclavable.
	D300 100 076	Laryngoscope 10mm 70° 159mm autoclavable.
	D300 100 077	Laryngoscope 10mm 90° 159mm autoclavable.
Otosopes		
	D300 100 078	Otoscope 1.9mm 0° 100mm autoclavable.
	D300 100 079	Otoscope 1.9mm 30° 100mm autoclavable.
	D300 100 080	Otoscope 2.7mm 0° 75mm autoclavable.
	D300 100 081	Otoscope 2.7mm 30° 75mm autoclavable.
	D300 100 082	Otoscope 2.7mm 70° 75mm autoclavable.
	D300 100 083	Otoscope 4mm 0° 50mm non-autoclavable.
	D300 100 084	Otoscope 4mm 30° 50mm non-autoclavable.

Photo	Reference	Description
Arthroscopes		
	D300 100 085	Arthroscope 1.9mm 0° 60mm autoclavable.
	D300 100 086	Arthroscope 1.9mm 30° 60mm autoclavable.
	D300 100 052	Arthroscope 2.7mm 0° 140mm autoclavable.
	D300 100 053	Arthroscope 2.7mm 30° 140mm autoclavable.
	D300 100 054	Arthroscope 2.7mm 70° 140mm autoclavable.
	D300 100 062	HD arthroscope 4mm 0° 185mm autoclavable.
	D300 100 063	HD arthroscope 4mm 30° 185mm autoclavable.
	D300 100 065	HD arthroscope 4mm 70° 185mm autoclavable.

Sheaths		
Photo	Reference	Description
For 2,9mm 12° or 30° hysteroscope		
	D300 110 002	Diagnostic sheath 4mm, rotatable stopcock.
	D300 110 000	Outer operating sheath 5.5mm, rotatable stopcock. To associate with D300 110 001.
	D300 110 001	Inner operating sheath 4.9mm, 7Fr. operative canal. To associate with D300 110 000.
	D300 110 032	Outer operating sheath, 5mm, rotatable stopcock. To associate with D300 110 034.
	D300 110 034	Inner operating sheath 4,3mm, 5Fr. operative canal. To associate with D300 110 032.
	D300 130 085	Examination bridge to associate with outer sheath of 18,5Fr. resectoscope, and with 5Fr. operative canal.


Instrumentation			
Photo	Reference	Description	
With flushing port			
 	D300 120 008	7Fr.	Graduated cissors, semi-rigid, L: 340mm. Single action jaw.
	D300 120 012	5Fr.	
 	D300 120 009	7Fr.	Graduated forceps (alligator), semi-rigid, L: 340mm. Double action jaws
	D300 120 013	5Fr.	
 	D300 120 010	7Fr.	Graduated biopsy forceps, semi-rigid, L: 340mm. Double action jaws
	D300 120 014	5Fr.	
 	D300 120 011	7Fr.	Graduated IUD grasping forceps, semi-rigid, L: 340mm. Double action jaws.
	D300 120 015	5Fr.	
Without flushing port			
 	D300 120 000	7Fr.	Graduated cissors, semi-rigid, L: 340mm. Single action jaw.
	D300 120 004	5Fr.	
 	D300 120 001	7Fr.	Graduated forceps (alligator), semi-rigid, L: 340mm. Double action jaws
	D300 120 005	5Fr.	
 	D300 120 002	7Fr.	Graduated biopsy forceps, semi-rigid, L: 340mm. Double action jaws
	D300 120 006	5Fr.	
 	D300 120 003	7Fr.	Graduated IUD grasping forceps, semi-rigid, L: 340mm. Double action jaws.
	D300 120 007	5Fr.	

Resectoscopes		
Photo	Reference	Description
18.5Fr. mini resectoscope		
	D300 130 000	Outer sheath, 18.5Fr., double flow.
	D300 130 001	Inner sheath, rotatable, 17.5Fr.
	D300 130 002	Passive working element for 18.5Fr. resectoscope.
	D300 130 009	Active working element for 18.5Fr. resectoscope.
	D300 130 003	Standard obturator for 18.5Fr. resectoscope.
26Fr. resectoscope		
	D300 130 004	Outer sheath, 26Fr., double flow.
	D300 130 005	Inner sheath, rotatable, 24Fr.
	D300 130 006	Passive working element for 26Fr. resectoscope.
	D300 130 008	Active working element for 26Fr. resectoscope.
	D300 130 007	Standard obturator for 26Fr. resectoscope.
	D300 130 010	Visual obturator for 26Fr. resectoscope.
	D300 130 063	Quick-lock adapter for Ellik evacuator connection. Compatible with MaxiCare & ReseCare.

Electrodes		
Photo	Reference	Description
18.5Fr. mini resectoscope		
	D300 140 000	REUSABLE
	D300 140 067	STERILE, SINGLE USE (Only by 10)
		Cutting loop 30°, resectoscope 18.5Fr.
	D300 140 001	REUSABLE
	D300 140 068	STERILE, SINGLE USE (Only by 10)
		Cutting loop, straight, resectoscope 18.5Fr.
	D300 140 002	REUSABLE
	D300 140 069	STERILE, SINGLE USE (Only by 10)
		Electrode knife, resectoscope 18.5Fr.
	D300 140 003	REUSABLE
	D300 140 070	STERILE, SINGLE USE (Only by 10)
		Electrode rollerball, resectoscope 18.5Fr.
26Fr. resectoscope		
	D300 140 004	Cutting loop 30°, resectoscope 26Fr.
	D300 140 005	Electrode knife, resectoscope 26Fr.
	D300 140 006	Electrode rollerball, resectoscope 26Fr.
	D300 140 007	Cutting loop, straight, resectoscope 26Fr.
	D300 140 008	Vaporisation electrode - rollerball with tip, resectoscope 26Fr.
	D300 140 009	Vaporisation electrode - rollerball smooth, resectoscope 26Fr.
Operative sheath		
	D300 140 011	Bipolar electrode, flexible, needle-type, 5Fr.

Photo	Reference	Description
	D300 140 010	Bipolar electrode, flexible, button-type, 5Fr.
	D300 140 038	Bipolar electrode, flexible, curved needle-type, 5Fr.

Light cables		
Photo	Reference	Description
	D200 150 500	Light cable, grey, Ø: 3.5mm, L: 2.3m with Storz connectors.
	D200 150 000	Light cable, grey, Ø: 3.5mm, L: 2.3m.
	D200 150 001	Connector Storz for light cable - endoscope side.
	D200 150 002	Connector Storz for light cable - light source side.
	D200 150 003	Connector Olympus for light cable - light source side.
	D200 150 004	Connector ACMI / Stryker for light cable - light source side.

Sterilization		
Photo	Reference	Description
	D200 160 500	Sterilization basket 480x250x47mm with a silicone mat & a lid.
	D200 160 001	Blue silicone mat 440x230mm.


delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

Photo	Reference	Description
	D200 160 004	Sterilization basket 480x250x47mm.
	D200 160 005	Lid for sterilization basket, 481x254mm.
	D200 160 502	Sterilization basket for 2 instruments.
	D200 160 003	Sterilization tube for electrodes.


Miscellaneous		
Photo	Reference	Description
Bipolar resection cable (ReseCare & MaxiCare)		
	D200 190 000	Covidien bipolar resection cable. Compatible with the Valleylab ForceTriad generator.
	D200 190 001	Erbe bipolar resection cable. Compatible with Erbe Vio 300 and Storz Autocon III 400 generators.
	D200 190 008	Bowa bipolar resection cable. Compatible with Bowa Arc 350 and 400 and Storz Autocon III 400 generators.
	D200 190 006	2 pins - 29mm bipolar resection cable. Compatible with all 2 pins - 29mm connection generators such as Erbe ES350, Pulsar MB 350, Megadyne and Lamidey Noury.


Photo	Reference	Description
	D200 190 007	Emed (SDS plug) bipolar resection cable. Compatible with the Emed Spectrum generator.
	D200 190 012	Berchtold or Martin bipolar resection cable. Compatible with Berchtold, Aesculap and Martin Maxium generators.
Monopolar resection cable (ReseCare)		
	D200 190 005	Valleylab monopolar resection cable. Compatible with Bovie, Valleylab and Conmed generators.
	D200 190 004	Erbe T-series monopolar resection cable. Compatible with Erbe T-serie, Martin and Berchtold generators (banana plug).
	D200 190 011	Berchtold monopolar resection cable. Compatible with Martin, Berchtold, Aesculap and Rudolf generators.
	D200 190 010	Erbe ACC/ICC monopolar resection cable. Compatible with Erbe ACC/ICC and Storz generators.
HF cable for bipolar electrodes (EasyCare)		
	D200 190 003	Valleylab cable for bipolar electrode. Compatible with Valleylab generators and international models.
	D200 190 002	Erbe cable for bipolar electrode. Compatible with Erbe and Storz generators.
	D200 190 009	Berchtold cable for bipolar electrode. Compatible with Berchtold and Martin generators.

Photo	Reference	Description
Other		
	D300 195 000	Reusable sealing cap for operative sheath.
	D300 195 001	Single use sealing cap for operative sheath.
	D300 190 000	Sterilization tube for endoscopes.

Uterine manipulators		
Photo	Reference	Description
Endojack		
	D250 100 000	Box of 5 steriles ENDOJACK uterine manipulators. 60-70mm size.
	D250 100 001	Box of 5 steriles ENDOJACK uterine manipulators. 80-90mm size.
Mangeshikar		
	D300 120 100	Handle with Pozzi forceps rotation module for Mangeshikar uterine manipulator
	D300 120 101	Outer sheath with fixation nut for Mangeshikar uterine manipulator
	D300 120 102	Inner sheath with locking pin for Mangeshikar uterine manipulator
	D300 120 103	Sliding sheath to use with cap for Mangeshikar uterine manipulator
	D300 120 107	Pozzi forceps for cervical holding, for Mangeshikar uterine manipulator

Photo	Reference	Description
	D300 120 104	Cap of 35mm with white insulating ring
	D300 120 105	Cap of 40mm with white insulating ring
	D300 120 106	Cap of 45mm with white insulating ring
	D300 120 108	Intra-uterine insert 4x40mm
	D300 120 109	Intra-uterine insert 5x50mm
	D300 120 110	Intra-uterine insert 6x60mm
	D300 120 111	Intra-uterine insert 7x70mm
	D300 120 112	Intra-uterine insert 8x80mm
	D300 120 113	Intra-uterine insert 8x120mm
	D300 120 114	Intra-uterine insert 8x150mm

3. Advertising items






Leaflets and booklets		
Photo	Reference	Description
	D900 700 055	EasyCare leaflet. French version.
	D900 700 010	EasyCare leaflet. English version.
	D900 700 054	ReseCare leaflet. French version.
	D900 700 009	ReseCare leaflet. English version.
	D900 700 013	MaxiCare leaflet. French version.
	D900 700 020	MaxiCare leaflet. English version.
	D900 700 011	iCare+ leaflet. French version.
	D900 700 019	iCare+ leaflet. English version.

Photo	Reference	Description
	D900 700 062	iFlow leaflet. French version.
	D900 700 063	iFlow leaflet. English version.
	D900 700 028	Endojack leaflet. French version.
	D900 700 029	Endojack leaflet. English version.
	D900 700 032	EasyCare clinical booklet. French version.
	D900 700 033	EasyCare clinical booklet. English version.
	D900 700 065	ReseCare clinical booklet. French version.
	D900 700 066	ReseCare clinical booklet. English version.

Accessories		
Photo	Reference	Description
	D900 700 058	Delmont imaging surgical cap.
	D900 700 059	Delmont imaging USB flash drive, 4Gb.
	D900 700 057	Delmont imaging purple canvas bag.
	D900 700 056	Delmont imaging retractable ballpoint pen.