

Resection electrodes

Photo	Reference	Description
18.5Fr. mini resectoscope		
	D300 140 000	REUSABLE
	D300 140 067	STERILE, SINGLE USE (Only by 10)
	D300 140 001	REUSABLE
	D300 140 068	STERILE, SINGLE USE (Only by 10)
	D300 140 002	REUSABLE
	D300 140 069	STERILE, SINGLE USE (Only by 10)
	D300 140 003	REUSABLE
	D300 140 070	STERILE, SINGLE USE (Only by 10)
26Fr. resectoscope		
	D300 140 004	Cutting loop 30°, resectoscope 26Fr.
	D300 140 005	Electrode knife, resectoscope 26Fr.
	D300 140 006	Electrode rollerball, resectoscope 26Fr.
	D300 140 007	Cutting loop, straight, resectoscope 26Fr.
	D300 140 008	Vaporisation electrode - rollerball with tip, resectoscope 26Fr.
	D300 140 009	Vaporisation electrode - rollerball smooth, resectoscope 26Fr.

Note: minimum delivery quantity of 5 reusable electrodes or one box of 10 single use electrodes.