




imaging all the women



NEW

MaxiCare

The 9mm bipolar resection

MaxiCare references and features:

HD Hysteroscope, 4mm, 12°

Autoclavable, L: 300mm

Ref. D300 100 002



Outer sheath, 26Fr.

Double flow. Quick lock.

Ref. D300 130 004



Inner sheath, 24Fr.

Rotatable


Ref. D300 130 005



Passive working element

Compatible bipolar

Ref. D300 130 006



Light cable

Autoclavable, Ø: 3.5mm, L: 2.3m. Supplied with Storz connectors

Ref. D200 150 500



NEW

Electrodes:

Description	Reference	
	Sterile, single use (Delivery in boxes of ten)	Reusable
 Cutting loop, 30°	D300 140 078	D300 140 004
 Cutting loop, straight	D300 140 081	D300 140 007
 Coagulating electrode knife	D300 140 079	D300 140 005
 Coagulating electrode rollerball	D300 140 080	D300 140 006
 Vaporisation electrode rollerball with pikes	D300 140 082	D300 140 008

Accessories:

NEW

MaxiCare sterilization basket

Autoclavable. 480x240x70mm. With personalized silicone holders

Ref. D200 160 505



Sterilization tube

Autoclaving and storage of electrodes

Ref. D200 160 003



Obturator

Quick lock

Ref. D300 130 007

Visual obturator

Quick lock

Ref. D300 130 010

Covidien bipolar resection cable

ForceTriad

Ref. D200 190 000

Erbe bipolar resection cable

VIO300D

Ref. D200 190 001



Zone Athélia V - 390, Avenue du Mistral - 13600 La Ciotat - France
Tel : +33 (0) 9 51 51 30 30 - Fax : +33 (0) 9 57 51 31 00
www.delmont-imaging.com - contact@delmont-imaging.com

Delmont Imaging is the manufacturer of the optic, Class IIa medical device, and of the resectoscopy system of Class IIa/IIb, marked CE 1304, notified body SIO, presented in this document. Delmont Imaging is the manufacturer of the light cable, Class I medical device.

Non contractual document - Réf. D900 700 020 D - 2020-03 - Copyright © 2017 DELMONT IMAGING. All rights reserved. No information or part of this document may be reproduced or transmitted in any form without the prior permission of DELMONT IMAGING