

Datasheet

Single use HF electrodes

General informations	
Conformity	Conform to 93/42/CEE
Applicable standards	See CE conformity Declaration for more details
Toxic/allergenic products	Latex free, DEHP free, Halogen free

HF electrodes			
Regulation			
Description	HF electrodes for resection system		
Intended use	Gynecological endoscopy only. Tissue removal, cutting, coagulation.		
Photos	 Cutting loop, 30° Cutting loop, straight Electrode rollerball Electrode knife Vaporisation electrode, rollerball with tip		
Device class	IIb		
Single-use electrodes			
Name	Cutting loop, 30°	Cutting loop, straight	Electrode knife
Article reference	D300 140 067	D300 140 068	D300 140 069
GMDN code	61872		
GTIN/UID ID code	3701217802766	3701217802773	3701217802780
Name	Electrode rollerball	Vaporisation electrode, rollerball with tip	
Article reference	D300 140 070	D300 140 091	
GMDN code	61872		
GTIN/UID ID code	3701217802797	3701217807051	

Technical data	
Colour code	Purple
Compatibility	18,5Fr. resectoscope
Generator output power	max. 2 kVp
Recommended power cutting mode	120-180 W
Recommended power coagulation mode	max. 100 W
Monopolar solution to be used	Glycine or Purisole for example
Bipolar solution to be used	0.9% NaCl
Packaging	10 units / box
Autoclavable	None. Already sterile and single use
Sterilisation method	ETO

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411