

Company							
Product name	NeoCare	Bigatti		Omni		Piranha	Mini Operating Hysteroscope
Photo	See annex of this document						
Hysteroscope							
Diameter	2.8mm	NP		NP		NP	2mm
Optic system	lenses	optic fiber		lenses		lenses	optic fiber
Angulation (DOV)	12°	6°		0°		0° or 12°	12°
FOV	80°	90°		80°		NP	90°
Diagnostic sheath							
Diagnostic sheath	✗	✗		✓		✗	✗
Outer diameter	/	/		3,7mm		/	/
Outer operative sheath							
Outer diameter	9mm	6.3mm	8mm	5.5mm	6mm	9mm	5.8mm
Operative canal diameter	4.5mm (13,5Fr.)	4mm (12Fr.)	NP	3mm (9Fr.)	4mm (12Fr.)	5mm (15Fr.)	3mm (9Fr.)
Working lenght	187mm	200mm	NP	200mm		NP	NP

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

*NP : Information not provided by the competitor

Company						
Product name	NeoCare	Bigatti		Omni	Piranha	Mini Operating Hysteroscope
Obturator						
Optional handle						
Dissociable sheaths						
Luer-Lock for cleaning						
Quick-Lock		monobloc				monobloc
Instrumentation						
Handle of the instruments						NP

delmont imaging


Athélia V - 390, Av. du Mistral - 13600 La Ciotat


Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00


www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

*NP : Information not provided by the competitor

Company						
Product name	NeoCare	Bigatti		Omni	Piranha	Mini Operating Hysteroscope
Diameter	3mm	4mm	3mm	/	/	3mm
Rigid instruments				/	/	
Cold instruments	 Straight micro-scissors	  	 Biopsy forceps	/	/	 Alligator forceps
	 Universal grasping forceps		 Alligator forceps			
	 Grasping and extracting forceps					

Company	 delmont imaging				
Product name	NeoCare	Bigatti	Omni	Pirahna	Mini Operating Hysteroscope
Morcellator		Shaver GYN  oval cutting window (4mm)  square cutting window (4mm)	Myosure procedure  Manual  Reach  XL (6mm sheath)		
Bipolar electrodes	 L-Hook bipolar electrode  Straight bipolar electrode  V-Hook bipolar electrode	 coagulation electrode			




HOLOGIC®
The Science of Sure



RICHARD
WOLF 

Piranha



Mini Operating Hysteroscope



delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

*NP : Information not provided by the competitor