

Competitive comparison

Hysteroscopy system

For hysteroscopes from 2.7 to 3 mm only.

Company		 SOPRO COMEG								
Product name	EasyCare	No name	No name	Bettocchi	OES 4000			Opti-Flow	No name	Gubbini Original
Photo	See annex of this document									
Hysteroscope										
Diameter	2,9 mm	2,7 mm	2,9 mm		3 mm			2,7 mm		2,9 mm
Angulation	30°	30°	12° and 30°	30°	30°			12° and 30°	25°	30°
Type of technology	Lenses	Lenses	Lenses		Lenses			Lenses		Lenses
Diagnostic sheath										
Diameter	4 mm	3,5 mm	4,1 mm		4 mm			3,5 mm		✗
Centimeter graduation	✓	✗	✗		✗			✗		
Outer operative sheath										
Diameter	5,5 mm	5 mm	5 mm	5,5 mm	4,5 mm	5,5 mm	6,5 mm	5,5 mm	5 mm	5,4 mm
Working length	216 mm	220 mm	210 mm	170 mm	187 mm			219 mm	205 mm	194 mm
Centimeter graduation	✓	✗	✓	✗	✗			✗		✗

For hysteroscopes from 2.7 to 3 mm only.

Company								
Product name	EasyCare	No name	No name	Bettocchi	OES 4000	Opti-Flow	No name	Gubbini Original
Atraumatic outer sheath	✓	✓	✗	✓	✗	✗	✓	✗
Flow	Double	Double	Double		Double	Double		Double
Dissociable outer and inner sheaths	✓	✓	✗	✓	✗	✗	✓	✓
Inner operative sheath								
Diameter	4,9 mm	NP	✗ mono-bloc	4,3 mm	4,9 mm	NP		NP
Centimeter graduation	✓	✗		✗		✗		✗
Atraumatic inner sheath	✓	✗		✗		✗		✗
Operative canal diameter	7 Fr.	5 Fr.	5 Fr.		7 Fr.	3 Fr.	5 Fr.	7 Fr.
Length of the insertion of the operative canal	Long	Short	Short	Long		Short		Long
								Short

For hysteroscopes from 2.7 to 3 mm only.

Company										
Product name	EasyCare	No name	No name	Bettocchi	OES 4000			Opti-Flow	No name	Gubbini Original
Rotatable stopcocks	✓	✓		✗	✗				✗	✗
Quick-lock between the two sheaths	✓	✗	✗ mono-bloc	✗	✗			✗ mono-bloc	✗	✓
Compatibility with other endoscopes (non captive connection)	✓	✓	✓		✗			✓		✓
Bipolar electrode compatibility	✓	✓	✓		✗	✓		✓		✓
Instrumentation										
Number of instruments proposed	4	5	10	3	1	5		4	3	6
Diameter	7 Fr.	5 Fr.	5 Fr.	7 Fr.	3 Fr.	5 Fr.	7 Fr.	5 Fr.	7 Fr.	5 Fr.
Working length	340 mm	400 mm	340 et 400 mm	400 mm	NP			340 mm	355 mm	335 mm










For hysteroscopes from 2.7 to 3 mm only.

Company								
Product name	EasyCare	No name	No name	Bettocchi	OES 4000	Opti-Flow	No name	Gubbini Original
Centimeter graduation					NP			
Rotation module					NP			
Luer-Lock for cleaning	 as desired				NP			
Quick-Lock					NP			
Handle of the instruments					NP			

Competitive comparison

Hysteroscopy system


For hysteroscopes from 2.7 to 3 mm only.

Company		 SOPRO COMEG							
Product name	EasyCare	No name	No name	Bettocchi	OES 4000	OptiFlow	No name	Gubbini Original	
Dismantable handle					NP				
Instruments	 Scissors  Biopsy forceps  Forceps  Grasping forceps	 Scissors  Biopsy forceps  Forceps  Grasping forceps	 Scissors  Scissors  Biopsy forceps  Biopsy forceps	 Scissors  Biopsy forceps  Forceps	Photos not provided		 Scissors  Biopsy forceps  Biopsy forceps  Forceps	 Hooked scissors  Biopsy forceps  Biopsy and grasping forceps	 Scissors  Biopsy forceps  Biopsy forceps
					Biopsy forceps	Scissors Hooked scissors Biopsy forceps Forceps with smooth and pointed jaws			

Competitive comparison

Hysteroscopy system

For hysteroscopes from 2.7 to 3 mm only.

Company		 SOPRO COMEG	 KARL STORZ ENDOSKOPE			
Product name	EasyCare	No name	No name Bettocchi	OES 4000	Opti-Flow No name	Gubbini Original
Instruments			 Forceps  Punch forceps  Instrument for myoma fixation  Palpator			 Forceps  Grasping forceps  Instrument for myoma fixation

For hysteroscopes from 2.7 to 3 mm only.

Annex:

