

Company														
Product name	ReseCare	MaxiCare	No name	No name			TCRis		Princess		Resec- tion Master	Gubbini		
												Ellipse	Original	Midi
Photo	See annex of this document													
Size of the corresponding hysteroscope	2.9 mm	4 mm	4 mm	2.9 mm	4 mm	4 mm	2.7 mm	3.3 mm	2.9 mm					
Resectoscope diameter	18.5 Fr.	26 Fr.	26 Fr.	15 Fr.	22 Fr.	26 Fr.	25.5 Fr.	21 Fr.	24.5 Fr.	27 Fr.	15 Fr.	16 Fr.	18.5 Fr.	
External shape of the resectoscope	Round		Round	Round			Round	Round			Ellipse	Round		
Working element	Passive or active		Passive or active	Passive			Passive or active	Passive or active	Pas- sive	Passive	Pas- sive or active			
Dissociable sheaths	✓		✓	✓			✓	✓			✓			
Rotative sheaths	✓		✓	✓	As desired		As desired	As de- sired	✓	NP*	As desired			
Quick-lock between the two sheaths and the working element	✓		✓	✓			✓	✓			✓			


Company												
Product name	ReseCare	MaxiCare	No name	No name			TCRis	Princess	Resec- tion Master	Gubbini		
										Ellipse	Original	Midi
Compatibility with other endoscopes (non captive connection)	✓		✓	✓			✗	✓		✓		
Resectoscope and electrodes usable in monopolar	✓		✗	✗			✗	✗		✓		
Automatic chip suction	✗		✗	✗			✗	✗	✓ but works poorly	✗		
Isolation	Black ceramic		Black ceramic	White ceramic			Black ceramic	Black ceramic		Black ceramic		
Number of electrodes offered	5	11	5	1	3	5	9	4	3	3	5	5
Cutting loop electrode, angled	R*+SU*	R+SU	R	SU	R	R 2 sizes	R 6 sizes	SU	SU	R+SU	R+SU 2 sizes	R
Cutting loop electrode, straight	R+SU	R+SU	✗	✗		R	✗	✗		✗	R+SU	R
Knife electrode	R+SU	R+SU	R	✗	R	R	R	R	R	R+SU		R

*R : Reusable electrodes

*SU : Single use electrodes

Company			 SOPRO COMEG	 KARL STORZ ENDOSKOPE								
Product name	ReseCare	MaxiCare	No name	No name			TCRis	Princess	Resec- tion Master	Gubbini		
										Ellipse	Original	Midi
Rollerball electrode	R+SU	R+SU	R	✗			✗	SU	SU	✗		R
Ball electrode	✗	✓	R 2 sizes	✗	R	R	R	✗		R+SU		R
Vaporization electrode	R+SU	R+SU 6 types	✗	✗			R	SU	✗	✗		
Comments				Unusual shape, chips get stuck in the elec- trode. 	Return bar on the electro- des, the chips get stuck in the electrode. 			Only single use, ex- cept the knife.		Unusual shape, chips get stuck in the electrode.		
										    		
										0° resectoscope.		

Annex:




ReseCare




SOPRO COMEG



OLYMPUS



TCRis



Princess 21Fr.

TONTARRA



Gubbini Original