

Competitive comparison

Bariatric resectoscopy system

Company	 delmont imaging	 STORZ KARL STORZ ENDOSKOPE	 OLYMPUS
Product name	ReseCare XXL	Resectoscope extra long	Long OES Pro Resectoscope
Photo	See annex of this document		
Claimed specialty	Gynecology	Urology	Urology
Size of the corresponding hysteroscope	2,9 mm	3,5 mm	4 mm
Working element	285 mm	265 mm	262 mm
Resectoscope diameter	18,5 Fr.	26 Fr.	27 Fr.
External shape of the resectoscope	Round	Round	Round
Working element	Passive or active	Passive or active	Passive or active
Dissociable sheaths	✓	✓	✓
Rotative sheaths	✓	✓	✗
Quick-lock between the two sheaths and the working element	✓	✓	✓
Bipolar resectoscope	✓	✓	✗


delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

Company	 delmont imaging	 STORZ KARL STORZ ENDOSKOPE	 OLYMPUS
Product name	ReseCare XXL	Resectoscope extra long	Long OES Pro Resectoscope
Also compatible with monopolar energy	✓	✗ (two different products)	/
Isolation	Black ceramic	Ceramic	White ceramic
Number of electrodes offered	5	3	3
Cutting loop electrodes, angled	R	R+UU	R+UU
Cutting loop electrodes, straight	R	✗	R+UU
Knife electrode	R	R+UU	✗
Rollerball electrode	R	✗	R+UU
Ball electrode	✗	✓ (monopolar)	✗
Vaporization electrode	R	✓ (bipolar)	✗
Other available product	Diagnostic set (30° hysteroscope and diagnostic sheath)	Hysteroscope and diagnostic sheath	✗
Comments	Unique on the gynecological market	Return bar on the bipolar electrodes, the chips get stuck in the electrode. 	/

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat


Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

*R : Reusable electrodes
*UU : Single use electrodes

Annex:



ReseCare XXL



Resectoscope extra long



Long OES Pro Resectoscope

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411