

Data sheet NeoCare



1304

The complete NeoCare set contains:

- 1 HD operative hysteroscope
- 1 outer sheath
- 4 plastic handles
- 4 3mm inserts
- 3 3mm bipolar electrodes
- 1 light cable
- 1 sterilization basket (2 trays)

General information	
Conformity	Conform to MDD 93/42/CEE
Applicable standards	See CE Conformity Declaration for more details
Toxic/allergenic products	Latex free, DEHP free, Halogen free

Hysteroscope	
Regulation	
Description	HD 8 mm 12° operative hysteroscope
Intended use	Gynecological endoscopy only
Photo	
Device class	Ila
Article reference	D300 100 051
GMDN code	36628
GTIN/UID ID code	3701217805460
Technical data	
Resolution	HD

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

Hysteroscope diameter	8 mm
Angulation (DOV)	12°
Field of view (FOV)	80°
Type of technology	Lenses
Working lenght	190 mm
Operative canal's diameter	4,5 mm
Autoclavable	Yes. See instructions for use for further details
Flow	Simple (Irrigation stopcock)
Connection for light cable	ACMI type
Light cable adapter	Type Wolf, Storz, Olympus
Outer sheath	
Regulation	
Description	Hysteroscopic outer sheath
Intended use	Hysteroscopy
Photos	
Device class	Ilb
Article reference	D300 130 004
GMDN code	35301
GTIN/UID ID code	3701217801318
Technical data	
Outer sheath's diameter	26Fr.
Working length	187 mm
Quick lock	Yes
Flow	Double (suction and irrigation stopcocks)
Autoclavable	Yes. See instructions for use for further details
3mm cold instruments	
Regulation	
Description	Reusable surgical instruments: scissors and forceps
Intended use	Endoscopic surgery

Device class	IIb			
GMDN code	62061			
Handle				
Name	Plastic handle for 3mm instruments			
Photo				
Article reference	D300 120 504			
Material	Plastic			
Ratchet	No			
Rotation module	Yes			
GTIN/UID ID code	3701217805194			
Autoclavable	Yes. See instructions for use for further details			
3 mm inserts				
Name	Straight micro-scissors	Curved scissors, Metzenbaum type	Universal grasping forceps	Grasping and extracting forceps, 2/3 teeth
Photos				
Article references	D300 120 509	D300 120 507	D300 120 510	D300 120 515
GTIN/UID ID code	3701217805149	3701217805163	3701217805132	3701217805088
Autoclavable	Yes. See instructions for use for further details			
Technical data				
Action jaws	1	2	2	1
Cleaning luer-lock	Yes			
Flexibility	Rigid			
Recommended energy	Cold instruments			
Diameter	3mm			
Working lenght	330mm			
3mm bipolar electrodes				
Regulation				
Description	3mm bipolar electrodes			
Intended use	Gynecological endoscopy only			
Photos				

Device class	IIb		
GMDN code	61873		
Bipolar electrodes			
Name	L-Hook electrode	Straight electrode	V-Hook electrode
Article reference	D300 120 772	D300 120 773	D300 120 774
GTIN/UID ID code	3701217805415	3701217805408	3701217805392
Autoclavable	Yes. See instructions for use for further details		
Technical data			
Energy	Bipolar. Coagulation only		
Working lenght	360 mm		
Diameter	3 mm		
Flexibility	Rigid		
Reusable electrode	Yes		
Solution to be used	0.9% NaCl		
Light cable			
Regulation			
Description	Light cable		
Intended use	Endoscopy		
Photo			
Device class	I		
Article reference	D200 150 500		
GMDN code	35507		
GTIN/UID ID code	3701217800090		
Technical data			
Diameter	3,5 mm		
Lenght	2,3 m		
Colour	Grey		
Autoclavable	Yes. See instructions for use for further details		
Adapters supplied	Storz		

Sterilisation basket	
Regulation	
Description	2 trays sterilization basket
Intended use	Reprocessing
Photo	
Device class	I
Article reference	D200 160 515
GMDN code	13730
GTIN/UID ID code	3701217807891
Technical data	
Trays dimensions (WxHxD)	520x50x250mm
Tray weight without cover	1,7 kg
Cover weight	1,0 kg
Stackable trays and covers	Yes
Material	Trays: stainless steel Supports: black silicone
Autoclavable	Yes. <i>See instructions for use for further details</i>