


Data sheet

ReseCare XXL

The complete ReseCare XXL set contains:

- 1 HD bariatric hysteroscope
- 1 outer sheath
- 1 inner sheath
- 1 working element
- 5 HF bariatric electrodes
- 1 light cable



Passive set

Ref : D400 100 035

Code GTIN/UID ID : 3701217803725

Active set

Ref : D400 100 036

Code GTIN/UID ID : 3701217803732

Regulation	
Conformity	Conform to MDD 93/42/CEE
Applicable standards	See CE Conformity Declaration for more details
Toxic/allergenic products	Latex free, DEHP free, Halogen free

Hysteroscope	
Description	
Description	HD 2.9 mm 12° bariatric hysteroscope
Intended use	Gynecological endoscopy only
Photo	
Device class	Ila
Article reference	D300 100 087
GMDN code	36628
GTIN/UID ID code	3701217806474
Technical data	
Resolution	HD


delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411

Hysteroscope diameter	2.9 mm
Angulation (DOV)	12°
Type of technology	Lens
Working length	395 mm
Autoclavable	Yes. <i>See instructions for use for further details</i>
Connection for light cable	ACMI type
Light cable adapter	Wolf, Storz, Olympus type
Resectoscopy system	
Description	
Description	Set comprised of working element and associated sheaths
Intended use	Bipolar or monopolar resectoscopy
Photos	 Outer sheath Inner sheath Passive working element Active working element
Device class	Ilb
Outer sheath	
Article reference	D300 130 088
GMDN code	35301
GTIN/UID ID code	3701217806481
Outer sheath's diameter	18.5Fr.
Working length	287mm
Quick lock	Yes
Flow	Double
Autoclavable	Yes. <i>See instructions for use for further details</i>

Inner sheath		
Article reference	D300 130 089	
GMDN code	35301	
GTIN/UID ID code	3701217806573	
Inner sheath's diameter	17.5Fr.	
Isolation	Black ceramic (4.5mm)	
Rotatable sheath	Yes	
Separable internal and external sheaths	Yes	
Quick lock	Yes	
Autoclavable	Yes. <i>See instructions for use for further details</i>	
Working elements		
Action	Passive	Active
Article reference	D300 130 002	D300 130 009
GMDN code	37805	
GTIN/UID ID code	3701217801295	3701217801363
Compatible with bipolar energie	Yes	
Compatible with monopolar energie	Yes	
Quick-lock between the two sheaths and the working element	Yes	
HF electrodes		
Description		
Description	HF bariatric electrodes for resection system	
Intended use	Gynecological endoscopy only. Tissue removal, cutting, coagulation and vaporisation	
Photo	 Cutting loop 30° Cutting loop straight Rollerball electrode Knife electrode Vaporisation electrode - rollerball with tip	
Device class	IIb	

Reusable electrodes				
Name	Cutting loop 30°	Cutting loop straight	Knife electrode	Rollerball electrode
Article reference	D300 140 086	D300 140 087	D300 140 088	D300 140 089
GMDN code	61873			
GTIN/UID ID code	3701217806542	3701217806559	3701217806566	3701217806573
Name	Vaporisation electrode - rollerball with tip			
Article reference	D300 140 122			
GMDN code	61873			
GTIN/UID ID code	3701217803688			
Autoclavable	Yes. See instructions for use for further details			
Technical data				
Colour code	Purple			
Compatibility	Bariatric resectoscope 18.5Fr.			
Generator output power	max. 2 kVp			
Recommended power cutting mode	120-180 W			
Recommended power coagulation mode	max. 100 W			
Monopolar solution to be used	Glycine or Purisole for example			
Bipolar solution to be used	0.9% NaCl			
Light cable				
Description				
Description	Light cable			
Intended use	Endoscopy			
Photo				
Device classe	I			
Article reference	D200 150 500			
GMDN code	35507			
GTIN/UID ID code	3701217800090			

Technical data	
Diameter	3.5 mm
Lenght	2.3 m
Colour	Grey
Autoclavable	Yes. <i>See instructions for use for further details</i>
Connectors supplied	Storz

delmont imaging

Athélia V - 390, Av. du Mistral - 13600 La Ciotat

Tel. : +33(0) 9 51 51 30 30 - fax : +33(0) 9 57 51 31 00

www.delmont-imaging.com - contact@delmont-imaging.com

SAS au capital de 237 200 euros - RCS Marseille 817 710 411 - N° de TVA Intracommunautaire : FR 81 817 710 411