





2mm hysteroscope available!



As the care of all women has always been our first priority, a diagnostic set comes to complete our range of products.

This 2mm hysteroscope will provide you with extra serenity for your diagnostic hysteroscopies, even in challenging situations. Cervical stenosis, nulliparous, anxious or reluctant patients, reduce the risk of vasovagal episode after the examination.

Its reduced diameter facilitates the introduction and manipulation into the uterine cavity, while maintaining a high definition image quality. Diagnosis is easy and the hysteroscope provides a real comfort to use.

*This tool is the perfect complement to the 2.9mm hysteroscopes.
Together they will cover all the needs of your dedicated room!*

Compatible diagnostic and operative sheaths with this 2mm hysteroscope are also available. You will find below the references and descriptions of these products:

Reference	Description
D300 100 090	HD 2mm hysteroscope, autoclavable 260mm 30°
D300 110 025	Diagnostic sheath with rotatable stopcock, 3mm
D300 110 043	Inner operating sheath, 3.7mm, 5Fr. operative canal
D300 110 044	Outer operating sheath, 4.3mm
D300 120 012	Semi-rigid graduated scissors, 5Fr.
D300 120 013	Semi-rigid graduated forceps (alligator), 5Fr.
D300 120 014	Semi-rigid graduated biopsy forceps, 5Fr.
D300 120 015	Semi-rigid graduated IUD grasping forceps, 5Fr.