



imaging all the women



Gynko®

A new approach
to hysteroscopy


delmont
imaging

With **Gynko®**, keep it simple!

Several hysteroscopies with a single hysteroscope.

Gynko® is a sterile, single-use hysteroscopy sheath. It is completely sealed, which allows you to keep your endoscope sterile. Perform several diagnostic hysteroscopies without having to repeat the sterilization processes between each examination. Gynko® is also available with a camera cover.



An accurate picture.

The optical performances of the distal lens are similar to those of an endoscope. Keep your picture quality for an accurate diagnosis!

The ideal tool for the See&Treat.

Gynko® reduced diameter (only 3mm!) helps you performing diagnostic hysteroscopies.

The flexible working channel can be used for endometrial biopsy or removal of small polyps. The cervix is gradually dilated when inserting the instrument, for a real See&Treat application.

The operative sheath total diameter remains below 4mm even when using 7Fr. instruments!

An efficient surgical procedure that respects the patient!



References:

Gynko sterile sheath - 2,9mm 30° hysteroscope

7Fr. working channel - With camera plastic cover

Delivery by 10

Ref. D300 110 068

Gynko sterile sheath - 2,9mm 30° hysteroscope

7Fr. working channel - Without camera plastic cover

Delivery by 10

Ref. D300 110 069

Gynko sterile sheath - 2,9mm 30° hysteroscope

5Fr. working channel - With camera plastic cover

Delivery by 10

Ref. D300 110 062

Gynko sterile sheath - 2,9mm 30° hysteroscope

5Fr. working channel - Without camera plastic cover

Delivery by 10

Ref. D300 110 063