




imaging all the women



MaxiCare

The 9mm bipolar resection

MaxiCare references and features:

HD Cystoscope, 4mm, 12°

Autoclavable, L: 300mm

Ref. D300 100 012



Outer sheath, 26Fr.

Double flow. Quick lock.

Ref. D300 130 004



Inner sheath, 24Fr.

Rotatable

Ref. D300 130 005



Passive working element

Compatible bipolar

Ref. D300 130 006



Light cable

Autoclavable, Ø: 3.5mm, L: 2.3m. Supplied with Storz connectors

Ref. D200 150 500



Electrodes:

Description	Reference	
	Sterile, single use (Delivery in boxes of ten)	Reusable
 Cutting loop, 30°	D300 140 078	D300 140 004
 Cutting loop, straight	D300 140 081	D300 140 007
 Coagulating electrode knife	D300 140 079	D300 140 005
 Coagulating electrode rollerball	D300 140 080	D300 140 006
 Vaporisation electrode rollerball with pikes	D300 140 082	D300 140 008

Accessories:

MaxiCare sterilization basket

Autoclavable. 480x240x70mm. With personalized silicone holders

Ref. D200 160 505



Sterilization tube

Autoclaving and storage of electrodes

Ref. D200 160 003



Obturator

Quick lock

Ref. D300 130 007

Visual obturator

Quick lock

Ref. D300 130 010

Adapter

For Ellick evacuator connection

Ref. D300 130 063

Erbe bipolar resection cable

VIO300D

Ref. D200 190 001