

iFlow mini accessories, references and features:

iFlow mini. Irrigation pump. Hysteroscopy and Laparoscopy

Hysteroscopy program includes:
- Setting target pressure range: 60 to 200 mmHg
- Pressure measurement range: 60 to 250 mmHg ±15mmHg
Maximum suction power: 50kPa
Power supply: (100-230 V~ ; 50/60 Hz)
Power consumption: 36 VA
Dimensions: 185 x 230 x 130 mm
Weight: 2.0 kg
Leakage isolation : CF type

Complete set including:
- iFlow mini. Irrigation pump only
- 2 x Pressure cuffs (modele of your choice)
- Y-adapter for double pressure cuff connection
- Self-test tube for iFlow mini pressure sensor

Ref. D400 100 600



Irrigation pump only.
Ref. 110 100 012

Pressure cuffs:

3L cuff.
275 x 465 mm
Ref. D110 100 020

3L cuff. Long version.
275 x 535 mm
Ref



A constant pressure irrigation pump specifically designed for gynecology

Whether you use it in outpatient room or in operating theater, iFlow mini will adapt perfectly and meet all your needs. Its compactness makes it an ally of choice in tight spaces. In operating rooms, enjoy a complete irrigation solution, including the two essential programs in gynecology: hysteroscopy and laparoscopy. Remarkably easy to use, its key performances remain intact, making iFlow mini a major asset in all situations.



Technical expertise at the service of vision quality

iFlow mini is based on a simple and effective irrigation technology. The liquid bags are inserted into pressure cuffs. These are automatically inflated by the pump to ensure continuous flow and constant pressure, an essential component to obtain optimal image quality. Additionally, adjust patient height for a reliable pressure display and thus, obtain uterine dilation in line with your expectations.

